

## TABLE OF CONTENTS

CHAPTER 1 Introduction .....	5
1.1 Welcome .....	5
1.2 Limitation of Liability .....	5
1.3 Registration Card .....	5
1.4 Copyright .....	5
CHAPTER 2 Installing the new Program .....	7
2.1 Installing the program .....	7
2.2 Making Sure My computer Date is Setup Properly .....	7
2.3 Building a shortcut .....	7
2.4 Giving Users the Correct Access Rights .....	8
2.5 Uninstalling the Program .....	8
2.6 Updating from Earlier Versions .....	8
2.7 Network Installation .....	9
CHAPTER 3 - Quick Start .....	10
3.1 Setting up .....	10
3.2 The Speed Buttons .....	11
3.3 Moving from Field to Field .....	12
3.4 Entering a new Employee .....	12
3.5 Selecting Occupation, Department, Workgroup & Company .....	12
3.6 Conditional Entries .....	15
3.7 Entering a Test for that Employee .....	15
3.8 Test Type .....	15
3.9 Medical History .....	15
3.10 Entering the Thresholds Manually .....	15
3.11 Drawing the Thresholds with the Mouse .....	15
3.12 The Percentage Loss .....	15
3.13 Transferring Data from your Audiometer .....	16
3.14 Copying the Audiogram Graph to the Clipboard .....	16
3.15 Selecting a Tester .....	16
3.16 Selecting an Audiometer .....	16
3.17 Australian Standard Behaviour .....	18
3.18 Printing the Data .....	19
3.19 Adding a Scanned Tester Signature to the Audiogram Report .....	20
3.20 Printing Labels for One Employee .....	20
3.21 Printing a letter .....	20
CHAPTER 4 - Starting Testing Employees .....	24
4.1 Setting up Defaults to Simplify your Work .....	24
4.2 What Happens to the Defaults .....	24
4.3 A good habit .....	24
CHAPTER 5 - Things to Know .....	25
5.1 Selecting From a Combo Box .....	25
5.2 Selecting from a Tabbed Book .....	25
5.3 Selecting from a Date Field .....	26
5.4 Moving and Selecting in a Grid .....	26
5.5 Adding an Employee .....	26
5.6 Changing from Employee to Test .....	26
5.7 Adding a test to an Employee .....	26
CHAPTER 6 - The File Menu .....	29
6.1 Data Entry .....	29
6.2 Browse Employees .....	29

6.3 Browse Tests	29
6.4 Set User Defaults	29
6.5 Import from a File WA Text File	29
6.6 Run Another Program	29
6.7 Backup Data	30
6.8 Doing the Backup	30
6.9 Data File Statistics	31
CHAPTER 7 - Data Entry Reports	32
7.1 The Employee's Audiogram	32
7.2 The Form 411 Report (WA)	32
7.3 The CX Report	32
7.4 The Report Card	32
7.5 The Spirometry Report	32
7.6 The Percentage Loss Report	32
CHAPTER 8 - The Report Menu	35
8.1 Batch Audiograms	35
8.2 Batch Report Cards	35
8.3 Batch Form 411 (WA)	35
8.4 Batch CX Report	35
8.5 Victorian Regulation 2007 3.2.11 (Test shift)	35
8.6 Australian Standard Report	35
8.7 The Assessment Screen	36
8.8 Base Line Shift Report	36
8.9 Whole Person Impairment Report (Workcover NSW)	36
8.10 Whole Person Impairment Report (Comcare)	37
8.11 OSHA STS Report	37
8.12 Category Report	37
8.13 Design Your Own Reporting Criteria	37
8.14 Assessment Report	37
8.15 Commercial Vehicle and Train Driver Reports	37
8.16 Report Creation	37
8.17 Selecting Data for Reports	37
8.18 Employee Status Statistics	37
8.19 Age Statistics	37
8.20 Employee Medical Statistics	38
8.21 Test Parameters & Medical Statistics	38
8.22 Potential Liability	38
CHAPTER 9 - The List Menu	41
9.1 List Employees	41
9.2 Listing Output Format	41
9.3 Labels	41
9.4 Mail Merge	42
9.5 Creating a Mail Merge Letter with MS Word	42
9.6 Printing Mail Merge Envelopes with MS Word	43
9.7 List Tests	44
9.8 Export Data to a WA Text File	44
9.9 Export to Excel	44
9.10 Re-test Management	45
9.11 Re-test Scheduling by Departments in Company	45
9.12 List Companies to Contact For Retest	45
9.13 List Departments to Contact For Retest	45
9.14 Re-Train Management	46
9.15 Hearing Protection Management	46
CHAPTER 10 - The Utility Menu	47
10.1 Re-index File	47
10.2 Untag Employees	47

10.3 List Duplicates .....	47
10.4 List Deleted Records .....	47
10.5 Purge Deleted Records .....	47
10.6 Rename a Company / Department / Occupation .....	47
10.7 Update Department & Company Numbers .....	48
10.8 Update Employee Retest Date .....	48
10.9 Recalculate Percentage Loss .....	48
10.10 Switch a test from One Employee to Another .....	48
10.11 Restore Data .....	48
 CHAPTER 11 The Setup Menu .....	 49
11.1 Setting up a Password .....	49
11.2 Adding Users to the Security System .....	49
11.3 Entering Companies .....	50
11.4 Entering Departments .....	50
11.5 Entering Occupations .....	50
11.6 Entering Workgroups .....	50
11.7 Hearing Protection Devices .....	50
11.8 Setting up the Employee .....	53
11.9 Setting up the Testers .....	53
11.10 Setting up the Test .....	53
11.11 Setting up the Audiometer .....	54
11.12 Australian Standard .....	55
11.13 Potential Liability .....	55
11.14 Comments .....	55
11.15 Letters .....	56
11.16 Page Setup .....	56
11.17 Printer .....	57
11.18 Preferences .....	57
 CHAPTER 12 Setting Up the Audiometer .....	 58
12.1 General Procedure .....	58
12.2 My Laptop does not have a Serial Port .....	58
12.3 Installing a USB to Serial adapter .....	58
12.4 The OSCILLA SM910 / 930 / 950 and 960 .....	59
12.5 The Maico MA30 .....	59
12.6 The Maico MA15 .....	60
12.7 The InterAcoustics AD27 .....	60
12.8 The InterAcoustics AS 216 and AD 226 .....	60
12.9 The GSI 66 / 67 / 68 Audiometers .....	61
12.10 Audiometer Data Transfer .....	61
12.11 Reading from the Earscan or Microlab .....	62
12.12 Earscan / Microlab Cable Diagram .....	62
12.13 Using an Earscan with Air and Tymp .....	62
12.14 Using an Earscan or Microlab with Air only .....	63
12.15 Single test transfer with the Microlab .....	63
12.16 Multiple test transfer with the Microlab .....	63
12.17 Multiple Test Transfer with the MA800 / MA728 .....	64
12.18 Setting up the Microlab / MA728 / MA800 to store data automatically .....	64
12.19 Setting up the Maico MA800 .....	65
12.20 Setting up the Maico MA728m .....	65
12.21 The Rexton AS53, AS54 & DA63 .....	66
12.22 The Tracor RA 300 / 500 .....	66
12.23 The Tracor RA 400 .....	67
12.24 The AudioScan 1000 .....	67
12.25 The GN Resound Xeta & Danplex 300a .....	67
12.26 The AS 608, AD629e, Affinity, Calisto & Astera Audiometers .....	67
12.27 The MedRx Audiometers .....	68
12.28 The Aurical Audiometers .....	68
12.29 Trouble shooting Communication Problems .....	69

CHAPTER 13 The letter System .....	70
13.1 Principal of Operation .....	70
13.2 Creating / Altering the Form Letters .....	70
13.3 Inserting a Replaceable Field .....	70
13.4 Creating Gender Sensitive Sentences .....	70
13.5 Using Bold and Underline .....	73
13.6 Changing the Letter Margins .....	73
CHAPTER 14 The Spell Checker .....	74
14.1 Setting up the Spell Checker .....	74
14.2 Spell Checking a letter .....	74
CHAPTER 15 Principles of Operation .....	76
15.1 Percentage Loss .....	76
15.2 Presbycusis Correction .....	76
15.3 The % Loss Report .....	77
15.4 Normals Assessment .....	77
CHAPTER 16 Interpreting the Audio Test Results .....	78
16.1 Assessment of the Audiogram .....	78
16.2 Assessment of the Report Card .....	79
CHAPTER 17 - Network Use .....	81
17.1 Installation .....	81
17.2 Host Setup .....	81
17.3 Other computers Setup .....	81
17.4 Parameter Change .....	81
17.5 Exclusive Processes .....	82
17.6 Concurrent Processes .....	82
CHAPTER 18 - Interfacing with Other Softwares .....	83
18.1 Medical Director .....	83
18.2 MS Word .....	83
18.3 MS Excel .....	83
18.4 SAP .....	83
CHAPTER 19 The To-Do List .....	85
19.1 Adding a To-do .....	85
19.2 To-do Display .....	85
19.3 To-do Categories .....	85
19.4 Listing the To-dos .....	85
CHAPTER 20 The Security System .....	87
20.1 Enabling the Security System .....	87
CHAPTER 21 Headset Correction .....	89
21.1 What is Headset Correction .....	89
1.1 Enabling Headset Correction .....	90
APPENDIX A - Layout of Export/Import file .....	91
APPENDIX B - Layout of Update file .....	93
INDEX .....	95

# CHAPTER 1 Introduction

## 1.1 Welcome

---

Welcome to the **Workplace Audiometrics** software.

Included in your package you should have

- The Install CD
- The Start up Manual (you are reading it now)
- A Registration card

## 1.2 Limitation of Liability

---

Every effort has been made to make this Software and Manual as accurate and functional as possible. The liability of **Biotronic** will be limited to making available to purchasers, updates and corrections, as may be found necessary.

**Biotronic** reserves the right to revise this software and manual at any time without notice. In no event shall **Biotronic** be liable to you for any indirect, incidental, consequential, special or exemplary damages or data loss arising out of or in connection with your use or inability to use the Workplace Audiometrics software.

You must assume full responsibility for the selection of the Workplace Audiometrics software to achieve your intended purposes.

## 1.3 Registration Card

---

Please fill out the required information on your Registration card. Mail the upper part to **Biotronic** and keep the middle part for your records. The lower part can be returned with payment if you wish to extend the update and technical support for a year. Our support includes:

- Updates when legislation changes or new features are implemented for a period of six months after purchase or yearly by renewal.
- Receive news of updates and enhancements.
- Access to our technical help line.

In all communications please quote your Licence number (printed on the Registration card).

## 1.4 Copyright

---

This program is licensed for use on one single computer or one Local Area Network (LAN) for the network version. Multiple version or licences must be purchased for use on a Wide Area Network. This program is protected by copyright laws. It is against the law to make copies of this program and give, sell or otherwise transfer those copies to friends, associates or any other persons.

Sure, it's a great program, and they should have their own copies. But send them to Biotronic to buy their own copy. That's the right way to do it.





## CHAPTER 2 Installing the new Program

The **Workplace Audiometrics** software will work under Windows 95, 98, Windows 2000, Windows XP, Windows Vista and Windows 8 or 8.1 as well as inside Virtual servers.

The program requires a minimum screen Desktop area of 600 x 800. All computers today can display this resolution.

### 2.1 Installing the program

---

- Make sure all desktop programs are closed
- Insert the CD in the CD Rom drive and the installer will proceed automatically
- If no then click on "START", click on "RUN"
- In the Run window, type D:\SETUP, (where D: is your CD drive) then press Enter

The Setup program will run. You need to install the program in the suggested directory C:\AUDIO. The program will not work if you install it in the Program Files folder in Windows 7 and Windows 8. As well, this will make trouble shooting and future updates easier.

- The Setup will offer to insert a shortcut in the Start menu, click on Yes, so that you will be able to run the Workplace Audiometrics program (also called ScreenW) by clicking on Start and on the ScreenW Icon.

When the installation is completed you **must** make sure the correct date format is selected.

### 2.2 Making Sure My computer Date is Setup Properly

---

Selecting the country [English (Australia)] and the short date format [dd/MM/yyyy]  
Click on **Start , Settings, Control Panel**



#### **Windows 7**

Double click on **Region & Language** then on the **Format** tab



#### **Windows 8, 8.1 and 10**

Double click on **Region**

Format at the top should read English(Australia)

Short date MUST read: **dd/MM/yyyy**

If not, change to the above and click on **Apply**. You can now close the Control Panel.

#### **Important:**

All dates **MUST be entered as a 4 digit numbers**. Entering 99 instead of 1999 confuses the computer as there is no way for it to assess if you mean 1999 or 2099

### 2.3 Building a shortcut

---

The installation process will create a shortcut for you, but just in case:

- Right-click anywhere on the background of the Desktop
- Select New, then Shortcut

- either type "C:\AUDIO\SCREENW..EXE" or select "Browse" and select the SCREENW.EXE program from the C:\AUDIO folder where it resides.
- Click "next", select a Name for the Shortcut, type "ScreenW"
- Click "next"
- Make sure the working directory is the same as for the Cmd\_line in the PROGRAM tab.

Example:

Cmd\_line: C:\AUDIO\SCREENW.EXE  
Working : C:\AUDIO  
Batch file:  
Shortcut key: None  
Run : Normal Window

- Click on finish (or OK)

You can now run the **Workplace Audiometrics** program by clicking on Start and on the **ScreenW** icon.

**Important:** Never install the program in your C:\Program Files folder. It will not work.

---

## 2.4 Giving Users the Correct Access Rights

On some computers you will need to make sure the program has the permission to read and write to its files. If the program start as a Demo, then open Windows explorer, right click on Properties, select the Security tab and tick Full Control for the users.

---

## 2.5 Uninstalling the Program

The default installation directory is C:\AUDIO. All files in that directory can be deleted, the installer does not install any other files on your computer.

---

## 2.6 Updating from Earlier Versions

To upgrade from version 3.x or 4.x, the upgrade expects to find your old program and data in a directory called **C:\Audio**, this location can be changed, but if you are using the software on a network you will probably have to make a copy on a local drive in C:\Audio, then copy it all back.

The upgrade CD should start automatically otherwise, click on the Start button then Run, they type D:\Setup (where D: is your CD drive). The setup program will ask for the location of the program (C:\Audio is most likely). It will then proceed to install the new program. When this is done, it will automatically run the File Converter. The File Converter expects to find the previous version data files in the same directory.

When the File Converter has completed the conversion, you will have to close it manually by clicking on the cross at the top right of the File Converter window (light blue)

The Upgrade will install the new software in the same directory as the existing one, but will not alter any existing file. The existing files will be copied and converted to the new data format.



### **Running the new program**

Just click on the shortcut, or on **Start > Programs > ScreenW**





## Completing the Setup

The File Converter will attempt to convert all your settings, however depending on the previous version some settings will have to be re-entered.

To check this read the **Quick Start** section and the **Setup** section.

## 2.7 Network Installation

---

The software is only installed once. This program is used by all computers via the network

### Step 1:

Make sure that your network is setup. Select one computer on the network, for our purpose we will call it “The Server”, where the program will be installed. This computer must be accessible by other computer on the network, regardless if it is a real file server or not. The “Server” is the computer which will hosts the Workplace Audiometrics software and data.

### Step 2:

Install the software normally on the “Server”.  
When you insert the CD, the install program should run automatically, otherwise double click on the Setup.exe file.

**Step 3:** Share the folder where the program is installed. This is usually C:\Audio. To share a folder you need to use Windows Explorer, right click on the folder, select **Shares** and make sure that all users have **read and write** access to this folder. This is very important, otherwise the program cannot be used by other users.

### Step 4:

- Go to the computer or workstation you want to have access to the software.
- Using Windows Explorer, select **Tools > Map Network drive** or **Computer > Map Network drive**
  - In the Drive, select K: for example
  - In the Folder, select the folder on the server computer where you just installed the program
  - If the server does not appear in the list, you may need to contact your IT support person
  - make sure "Reconnect a logon" is ticked
  - Click on Finish
  - You can now access your server from this station using that drive letter

### Step 5:

- Creating a new Shortcut
- Right click the Desktop
  - Browse and locate the drive letter created above (K:)
  - Select the proper directory / folder
  - Click on the Software name (ScreenW.exe)

The Server computer must now always be turned ON first and turned OFF last, so that the software and data are always available to the other computers using the program.



## CHAPTER 3 - Quick Start

Run the program by clicking on the screen icon or click on **Start, Programs, ScreenW** or double clicking on the shortcut if you have built one.

We will now review how to configure the program setup to reflect some of your preferences. For a more complete overview, please refer to the **Setup** chapter.

Make sure that you click on **Save** after completing each window.

### 3.1 Setting up

---

- Go to the **Setup** menu, we will explore now the essential items required to get started. You should visit each item of this menu later to make sure the software is configured for your personal preferences.



#### **Your Company Details**

**Go to Setup > Preference > Company Details** to enter your details.

Depending on the country you select, country specific reports may become available relevant to the country you selected. There is also a Company title in **Setup > Hearing Test > Print out** which should be filled with your company name.



#### **Companies**

**Add** a new Company, or **Edit** an existing one to enter your own company details. Then insert that company name as the **Default Company**. If you work with more than one company, then you can enter each one. Click on **Exit**.



#### **Departments**

Add at least one Department, or modify an existing one. Departments can be added at the time they are selected in the Employee screen, when entering an Employee's data. So there is no need to enter all of them in advance. The number is updated automatically each time you select a Department. Click on **Exit**.



#### **Occupations**

As above.



#### **Workgroups**

A workgroup is an area of similar noise exposure in L Aeq, 8h. The noise level can be entered, so that when a test is added, the wearing of an adequate hearing protector can be checked for that person.

The software will automatically enter the corresponding noise class based on NAL and Australian Standard 1269, 1998. If the exposure is tonal (one specific frequency only) or above 105 dB this system is not applicable and a noise professional's help may be required. Click on **Exit**.



#### **Employee**

Click on the **Special Labels** tab and make sure the label names are suitable for your use. Each

label gives meaning to a field which contain your employee's data. So if the label says **Employee Number**, that is what is meant to appear in the field.

You may wish to change Employee Number to Bundy number for example. Click on **Save**.



### Test

In the **Test Setup** tab, select the **Percentage loss calculation formula** for your State. Since 2001 by law, It should be the 1988 tables except in Western Australia where the extended 1988 formula is used. Enter the Report titles to represent your company. Try a couple of printouts and alter the titles until you are happy with the result.

In the **Spiro / Tymp / Bone** tab, un-check the option you do not need. If you select the Bone conduction, you will be able to import Bone conduction data directly from an Oscilla 960 or AD27 audiometer. Click on **Save**.



### Testers

Go to **Setup > Testers** Click on **Add** and enter the name of the tester, add as many as required. Highlight the default tester and click on **Set Default**. Click on **Exit**.



### Audiometer

Go to **Setup > Audiometer**, Select the audiometer connected to your computer. If you requested it, the software is supplied with a cable for your audiometer providing it is supported by the software. Select the communication port, usually **Com1**. If you use a USB to Serial converter it may be **COM5** or even higher (I have seen Com20 being used).

Click on **Add** and enter the name, serial number and calibration date of all your audiometers. Highlight the default audiometer and click on **Set Default**. Click on **Exit**.



### Australian Standard

This feature allows you to store one test and one re-test done at the same time. You may not want to do tests this way, if so, turn the feature off. Click on **Save**.



### Preferences

If you wish to preview the reports before they print, Check **Preview Reports**, you can always turn it off later. Click on **Save**.

## 3.2 The Speed Buttons

---

There are 3 Speed Buttons.



### Employee

This button is replicated in the **File** menu. It is a quick way to browse through your employees. This option is not meant to be used for data entry, but data can be modified there. Employees are shown in Alphabetical or Work number order at your choice.

Employees can be tagged: click on Tag and select a tag. Tags can be customised in

**Setup, Employee**, click on the **status \$ Tags** tag

Employees have a status: to set your own, go to Setup, Employees, click on the **Status & Tags** tag.



### Tests

This button is replicated in the **File** menu. It is a quick way to browse through your employee's Test. This option is not meant to be used for data entry, but data can be modified there. Tests are shown in Entry Test Date order. The display can be set to only show All, Normal or Abnormal test and by a date range.

To set a date range: enter the "Date from". The "Date to" and click on **Search**

To cancel the date range: click on **Reset**.

To show only the normals: click on **Normal Tests only**

To cancel and show all again: click on **All Tests**



### Data entry

This button is replicated in the **File** menu. It is the quicker way to do all your data entry. Click on the **Data Entry** button on the main window. This is where all the data is entered.

---

## 3.3 Moving from Field to Field

The **TAB** key is used to move down to the next chronological field.

The **SHIFT TAB** is used to move back up. The mouse is an essential component to select fields at random. This could lead to not entering essential information. The software will therefore check that the essential information is entered and warn the user. In general fields highlighted in yellow are essential to the proper operation of the software.

---

## 3.4 Entering a new Employee

From the Main menu, Click on the **Data entry** button. Click on the **Employee Tab**. Click on the **Add button** (with a little person icon) to add a new Employee. The fields in yellow are essential to the workings of the program so they must be filled up. Remember to enter all dates with a 4 digit year eg 1999 or 2000, not 99 or 00, **this is very important**.

When the employee data entry is completed, click on **Save**. You can find an employee by name just by typing the name in the **Search** window (top right) the Employee grid will automatically find the name as you type.

You can change the display order from alphabetical to Employee number by selecting that order in the Order pull down box or by clicking on the column header of the column you want the order in. The corresponding column will be highlighted in light blue to reflect the order.

---

## 3.5 Selecting Occupation, Department, Workgroup & Company

- Click on the corresponding button
- To find an entry, click in the Search window and type the name, as you type the grid will zoom on the name in the list, when the name appears, just double click on the red highlight, or click on **Select**.

If you need to add a new entry, click on **Add**, enter the data, click on **Save**, then **Select**.



The main screen of the program

This is the screen you will see when starting the program. All data entry functions are performed by pressing the Data entry button.

The Employee and Test buttons allow browsing through lists of Employee and tests. The menu options allow access to all the reporting, listing, utility and setup functions.

Setup Employee screen layout and parameters

Screen Layout Employee Status & Training Custom Labels Import setup

**Design your own layout**

Title PREVIOUS EMPLOYERS

1 -  
2 -  
3 -

Title MEDICAL HISTORY

Years S

Exposure to explosive/gun fire  
Military service  
History of hearing troubles  
Exposure to childhood antibiotics

Noisy Hobbies  
Power Tools  
Dizziness  
Surgery  
Referred medically

Ringing ears  
Frequent Colds  
Discharges  
Head injury  
Mumps / Measles

\* To use the NAL 80 standard, This field must be called Discharge

Save Employee name in Uppercase  
Save Company name in Uppercase

**Final appearance**

PREVIOUS EMPLOYERS

1 -  
2 -  
3 -

MEDICAL HISTORY

Exposure to explosive/gun fire  
Military service  
History of hearing troubles  
Exposure to childhood antibiotics

Noisy Hobbies  
Power Tools  
Dizziness  
Surgery  
Referred medically

Ringing ears  
Frequent Colds  
Discharges  
Head injury  
Mumps / Measles

This panel shows what your design will look like.

Reset to Defaults Show Field names

See also the Test Setup page

Save Cancel

The Employee information setup screen

Setup test and audiogram display

Test setup Audiogram Report Card Test Types & Questions Spiro, Bone & Tym

**Percentage Loss Calculations**

Calculate Percentage Loss  
Correct % Loss for age

Calculation method 2- 1988 Tables (Standard ACT,NSW,VIC,QLD)

**Re-test**

These values are used to automatically set the re-test date, based on the date of test

Retest noise exposed Employees every 2 Years

Retest noise exposed Employees 0 months after baseline test (zero = ignore)

Retest high risk employee every 12 months

Retest Non-Exposed Employees every 2 Years (for example office workers)

**Printout**

Test Report Title Hearing Assessment Report

Test Report Subtitle Mobile Screening Pty Ltd

Percentage Loss Report Title Detailed Percentage Loss Report

Show company name on all reports

Name label for Baseline tests Baseline test

Name label for Routine tests Routine test

Save Cancel

The test information setup screen

### 3.6 Conditional Entries

---

Some entries are conditional to the settings of other fields. For example, you cannot enter a Termination date if the Employee status is not Terminated.

If **Allow Re-Test** or **Allow Re-Train** are checked, then **Next test due** and **Next training** will be calculated when you enter a test. The calculation will be based on the **Re-Test** and **Re-Train** gap you selected in the **Setup menu, Test, Test Setup** tab.

### 3.7 Entering a Test for that Employee

---

Click on the **Test** tab, click on **Add** (The button with the headphone icon)

### 3.8 Test Type

---

<p><b>Important:</b> The BASELINE Test is normally the first test chronologically. But it can be entered at any time so long the date is correct. The software will sort the test by date automatically.</p>
--

### 3.9 Medical History

---

The test contains the same history as the “Employee’s record.” This is to allow you to keep track of the initial history and of the changes occurring at the time of each test.

### 3.10 Entering the Thresholds Manually

---

You can transfer tests from your audiometer or enter the figures manually. Double click on the first threshold then type away, if the threshold has only one digit, enter a space afterwards, that way the cursor will jump automatically to the next threshold making data entry faster.

### 3.11 Drawing the Thresholds with the Mouse

---

Click on the Left or Right button, then click on the correct location on the graph. To erase a threshold, just click at the bottom below 95dB for that frequency. The threshold values will be automatically written in the correct location.

If you have selected to enter Bone conduction, you will have 2 extra buttons to draw Bone thresholds. No line will be drawn, for clarity, but the threshold values will be written in the correct locations on the Bone page.

This feature can be disabled in **Setup > Test > Audiogram** if necessary, in which case the Mouse Draw buttons will not be visible. You do not have to click on the **Off** button, it is only there to avoid changing thresholds by mistake.

### 3.12 The Percentage Loss

---

The percentage loss is calculated when you save the test, or when you click on the **Recalculate** button.

The Employee’s age is calculated from date of birth and date of the test. The length of employment is calculated from the date of commencement and the date of the test (not the date you enter the data, so that entering tests done at an earlier time will still have correct information).

When adding a test, the **Occupation** and the **Department** are transferred to the test automatically to keep a record of the Employee's environment at the time of the test.

### 3.13 Transferring Data from your Audiometer

---

If you have selected your audiometer as suggested in the Quick Setup, a **Capture Test** button will be visible in the test window. If you selected Australian Standard (or Esso Standard) a **Capture Re-Test** button will also be visible. Next to these buttons will be the name of your audiometer.

Just click on the button to access the Data Transfer window which will instruct you as to the specific procedure for transfer for this audiometer. If you have selected Australian Standard (or Esso Standard), each button will automatically transfer the data to the relevant thresholds.

**Important:** To select an Audiometer, go to the **Setup** menu, **Audiometer** and select the proper audiometer, serial port etc....

### 3.14 Copying the Audiogram Graph to the Clipboard

---

You may want to include the graph displayed in the Employee data entry window in a Word or Excel document. This is done by using the Windows Clipboard. The graph is copied to the Clipboard by clicking on the green button below the mouse on the left of the graph. To paste the graph, move to the document, position the cursor where you want the graph and press **Ctrl V** or click on **Edit > Paste**.

### 3.15 Selecting a Tester

---

The name of the tester can be selected from a pre-entered list. You can define this list in the **Setup > Tester** menu.

When you enter a new test, the software will automatically load the default Testing Officer, you can change this selection by clicking on the button on the right of the Tester field. The Testing Officer and the Approval number will be automatically loaded.

### 3.16 Selecting an Audiometer

---

The name of the Audiometer can be selected from a pre-entered list. You can define this list in the **Setup > Audiometer** menu. You can select the audiometer which is connected to your computer and the audiometers you may (or have used) to do tests.

When you enter a new test, the software will automatically load the default Audiometer, you can change this selection by clicking on the button on the right of the Audiometer field. The Audiometer, Serial number and calibration date will be automatically loaded.



Surname	Name	Work No	Status	Shift	Cost centre	Re test	Workgroup	Department
ASPARTATE	Garry		N			18/04/2011	OFFICE WORK	INFORMATION S
BELLA	DONNA		T	red	123456789*	12/06/2003	FOUNDRY	FOUNDRY
BLOB	Fred		W			01/08/2008	COOK	FOOD
BLOGGS	Freddy		P			17/12/2019	COOK	FOUNDRY
CAN	Larry	12313	H	ASD	123	05/02/2017	CLERK	FOUNDRY
DALLILA	Alison		N					PRESS
DAVIDSON	Jenny		H					PRESS
FILLMORE	Peter		P					PRESS
PAN	Peter		P			25/12/2017		PRESS
STRADIVARIUS	Viola	12345	C	45	12345	18/02/2007	CLERK	TYPING POOL
WALT	Disney		P			21/06/2017		PRESS
ZORRA	Guy	12	P	12		01/08/2007	FOUNDRY	PRESS

Employee	Test	To do
<b>CAN Larry</b> ID: LAR 10 <div> <span>Add</span> <span>Edit</span> <span>Del / Undel</span> <span>Print</span> </div>		
<b>Status: Permanent - High risk</b>		
Street 12 norfolk street Suburb SURU HILLS Phone 123 3456 State NSW Postcode 2121 Gender Male Date of Birth 01/01/1961 Entered 05/02/2002		
<b>Employment details</b> Work No 12313 Commencement date 21/12/2001 Cost centre 123 Occupation Serviceman Shift ASD Department FOUNDRY Company ROBERTSON AND SON Address 56 PARRAMATTA ROAD, Workgroup CLERK Noise Exposure Class 0 Protector Class 1		
<b>Comments</b> Works in high noise area		
<b>Actions allowed</b> <input checked="" type="checkbox"/> Re-test <input checked="" type="checkbox"/> Re-train		
<b>Further actions</b> Next test due on 05/02/2017 Next training due 05/02/2017		
<b>Claims</b> <input checked="" type="checkbox"/> Claim submitted <input type="checkbox"/> Claim settled		
<b>PREVIOUS EMPLOYERS</b> 1- sinkwell ship builders 2- ted's bakery 3- fanfare Pty Ltd		
<b>MEDICAL HISTORY</b> <input type="checkbox"/> Medical history changed <input checked="" type="checkbox"/> Referred medically <input checked="" type="checkbox"/> Noisy Hobbies <input checked="" type="checkbox"/> Power Tools <input checked="" type="checkbox"/> Dizziness <input checked="" type="checkbox"/> Head injury <input checked="" type="checkbox"/> Surgery <input checked="" type="checkbox"/> Exposure to explosive/gun fire 1/ears <input checked="" type="checkbox"/> Military service <input checked="" type="checkbox"/> History of hearing troubles <input type="checkbox"/> Exposure to childhood antibiotics <input checked="" type="checkbox"/> Ringing ears <input checked="" type="checkbox"/> Frequent Colds <input checked="" type="checkbox"/> Discharges <input type="checkbox"/> Mumps / Measles		
<b>History comments</b> Need medical referral		

The Data entry window with the Employee tab visible

The bottom part of the window shows all the information for Larry Can, the employee highlighted in red within the top part of the window. The status of this employee is displayed in a yellow box. Larry is permanently employed (not a contractor) and is exposed to high risk noise during the course of his work. Larry is tagged because he had a hearing loss claim. Tags can be configured to say anything you want to draw attention to.

Change an Employee	
Surname <b>CAN</b> *	Name <b>Larry</b> *
Work No 12313	
<div> <span>Save</span> <span>Cancel</span> </div>	
Street 12 norfolk street Suburb SURU HILLS Postcode 2121 State NSW Gender MALE * Telephone 123 3456	
<b>Actions allowed</b> <input checked="" type="checkbox"/> Allow Re-Test <input checked="" type="checkbox"/> Allow Re-train	
<b>Actions</b> Date of birth 01/01/1961 * Entered 05/02/2002 Next test due 05/02/2017 Next training 05/02/2017	
Occupation Select Occupation serviceman Department Select Department FOUNDRY Workgroup: CLERK Protector Class required 0 Company Select Company ROBERTSON AND SON 56 PARRAMATTA ROAD,	
<b>Comment</b> Works in high noise area <span>Spell check</span>	
<b>Workplace</b> Status Permanent - High risk * Cost centre: 123 Shift: ASD Commencement date 21/12/2001 Termination date <input checked="" type="checkbox"/> Claim submitted <input type="checkbox"/> Claim settled	
<b>PREVIOUS EMPLOYERS</b> 1- sinkwell ship builders 2- ted's bakery 3- fanfare Pty Ltd	
<b>MEDICAL HISTORY</b> <input type="checkbox"/> Medical history changed <input checked="" type="checkbox"/> Referred medically <input checked="" type="checkbox"/> Noisy Hobbies <input checked="" type="checkbox"/> Power Tools <input checked="" type="checkbox"/> Dizziness <input checked="" type="checkbox"/> Head injury <input checked="" type="checkbox"/> Surgery <input checked="" type="checkbox"/> Exposure to explosive/gun fire <input checked="" type="checkbox"/> Military service 1 Years <input checked="" type="checkbox"/> History of hearing troubles <input type="checkbox"/> Exposure to childhood antibiotics <input checked="" type="checkbox"/> Ringing ears <input checked="" type="checkbox"/> Frequent Colds <input checked="" type="checkbox"/> Discharges <input type="checkbox"/> Mumps / Measles	
<b>Comments</b> Need medical referral	


The fields in yellow are compulsory as they are required for calculations and assessments. All other fields are optional.

Actions allowed (re-test and re-train) control who is included in **Re-test** and **Re-train** management reports. These reports will provide information on who need to be re-tested or re-trained.

The name of the fields “**Emp No**”, “**Cost Centre**” and “**Shift**” can be altered to suit the need of your organization. They can be changed in **Setup > Employee > Custom labels**

### 3.17 Australian Standard Behaviour


When the Australian Standard is disabled the test results look like the image below: and the percentage loss is calculated using the thresholds displayed. All the report also use these thresholds.



When the Australian Standard is enabled The percentage loss is calculated using the thresholds in yellow in the image below. The following reports will also use these thresholds:

- Percentage loss report, - CX Report, - Audiogram report, - Report Card Report, - Graphic Report

All other reports use the results from the first test (the first two lines)



You can see that the percentage loss is slightly different due to the lower thresholds in yellow at 3K and 6K.

### 3.18 Printing the Data

---

All the reports will print to your default printer. They may display as Preview depending on the Preview flag, if so, click on the print button to send then to the printer.

To Preview a report, click on Preview before clicking on the relevant Print button. A setting is also available in the Setup menu for reports without a specific Preview button. To set this: go to **Setup**, select **Preference**, in the Report box, un-check **Pre-View Reports before Print**.



#### Printing the Audiogram

1- From the **Client** Tab

The latest Audiogram will automatically be printed when clicking on the **Print** button. A label indicates “Latest Test”.

2- From the **Test** Tab

The highlighted Audiogram will automatically be printed when clicking on the **Print** button. A label indicates “This Test”. In this way you can select which test you want to print simply by highlighting a different test.

Click on Preview to view the Audiogram on the screen first. You will then have an option to print it.



#### Printing the Percentage Loss

as for the Audiogram above, but select the % Loss button instead of the Audiogram button.



#### Printing the Audiometry Report Card

The Report card shows the Employee’s details and all the Audiometry tests available in the database for this Employee. The overall loss between the Base line (first test) and the last test is shown at the bottom of the report for each frequency. In the Data Entry window, click on the **Report Card** button, in the top right section( it does not matter if the Employee or the test tab is visible).



#### Printing the Spirometry Report Card

This option will not be enabled if you have not selected Spirometry in the **Setup > Test > Spiro/Tymp/Bone** menu option. The Spiro Report shows the Employee’s details and all the Spirometry tests available in the database for this Employee.

In the Data Entry window, click on the **Spiro Report** button , in the top right section( it does not matter if the Employee or the test tab is visible).



#### The assessment Report

This option will analyse the current Employee’s tests and report on 10 different criteria. 5 from the Australian Standards, 1 from the USA Rio Tinto Standard and 2 from the Victorian Regulation.

This report can be printed or exported. Two formats are supported: RTF (to open with MS Word) and Text suitable for import into SAP for example.



### **The Form 411 Report (for WA)**

This option will print a complete Form 411 for WA Workcover with all the details filled and an the employee's audiogram.



### **The CX Report**

A very good report to give to the person being tested. It shows the results with explanations in layman's terms.

---

## **3.19 Adding a Scanned Tester Signature to the Audiogram Report**

You can scan your signature and have it printed automatically on the Audiogram report. To do this, scan the signature and save it as a bitmap file (bmp) so that the length is 4 times the height. Save it with the same name as the tester and the software will automatically pick the right signature depending on the tester name.

Example: Fred Bloggs.bmp

---

## **3.20 Printing Labels for One Employee**

It is sometimes convenient to have the ability to print one or more labels for one Employee. For example, to use for an Employee report folder or to mail a letter. This feature will allow you to print any number of labels for the same Employee (the one highlighted in red in the Employee grid.)

You can select the label row and column to start printing so that the same sheet of labels can be re-used.

In the Data Entry window, click on **Print, Labels**. The label format can be selected on the **Setup, Page Format, Labels**. A number of preset formats are available or you can make up your own. Continuous labels are not supported. The first line of the labels can be set in Bold and in a larger font if required, for example for a folder label.

To print a label for each of a number of Employees, go to the **List** menu, and select **List Employee**, then make a group selection then click on **Label**, then **Start**. See the **List menu** for more information.

---

## **3.21 Printing a letter**

The software has a system of letters. The letters are created in the **Setup, Letters** menu. You can create as many letters as you wish. They are saved in RTF format(Rich Text) and can also be edited by MS Word.

The software supports a simple but powerful form of Mail Merge. Fields can be selected when the letter is created. They are replaced with the relevant employee data at the time of printing the letter. See the Chapter on Setting up letters for more information.



### **How Can I Find the Last Letter I Printed**

When opening a letter in the Data Entry, for example ABNORMAL.RTF, it is scanned for Mail Merge fields, the fields are decoded and the final product is shown to you so that you can alter it.

Employee and Test data

Surname	Name	Work No	Status	Shift	Cost centre	Re test	Workgroup	Department
ASPARTATE	Garry		N			18/04/2011	OFFICE WORK	INFORMATION S
BLOB	Fred		W			01/08/2008	COOK	FOOD
BLOGGS	Freddy		P			17/12/2019	COOK	FOUNDRY
CAN	Larry	12313	H	ASD	123	05/02/2017	CLERK	FOUNDRY
DALLILA	Alison		N					PRESS
DAVIDSON	Jenny		H					PRESS
FILLMORE	Peter		P					PRESS
PAN	Peter		P			25/12/2017		PRESS
STRADIVARIUS	Viola	12345	C	45	12345	18/02/2007	CLERK	TYPING POOL
WALT	Disney		P			21/06/2017		PRESS
ZORRA	Guy	12	P	12		01/08/2007	FOUNDRY	PRESS

Remember this Employee

Surname

Name

Order Employee by

Hearing loss claim ☐

Employee exposure Status

Permanent - Quiet environ

Terminated ☐

Copy

Assessment

Labels

Letters

Tag

Exit

Employee Test To do

Tested 14/12/2017 Routine test

Actions due

Next test 18/04/2011

Next training 18/04/2013

Hearing protection

Type of protection worn QB2

Hearing protection Class 2

☒ Protection worn before test

☒ Protection usually worn

☒ Training provided

Details at time of test

Employment 9 Years Age 62

Occupation ENGINEER

Department INFORMATION SYSTEMS

Audiometer Audiometer 1 Cal 3/04/2000

Tester Goody Gumdrop

Comments

Hearing Test Medical Spiro

Air Conduction thresholds

	.5k	1k	1.5k	2k	3k	4k	6k	8k	Hz
Left	15	25	25	25	25	45	55	50	dB
Right	10	15	15	15	10	15	20	30	dB

Audiogram

This test is **Abnormal** Left **Normal** Right

Visual examination normal normal

Driver test 22.5 dB 12.5 dB

% Loss 6.8 % 0.0 %

% Loss Binaural 1.6 %

% Loss age corrected, NAL Tables 1988

Headset correction ☐

The data entry screen showing the employee page

Print

➔

Audiogram

Form 411 (WA)

CX Report

Report card

Analysis

Spiro Report

% Loss

Output

☐ Print

☒ Preview

☐ pdf file

☐ rtf file

☒ Air ☐ Bone

☒ Standard ☐ Insurance

Select

Close

Printer=FS-41000N

The Print options window is activated when pressing the Print button. It give access to 4 reports which can be sent to the screen, the printer or created as a pdf file to attach to emails and share.

When the Test tab is visible, the Print button allow printing of any tests. When in the Client tab, only the latest test is printed, because it is not possible then to select which test to print

Spirometry, Tympanometry and Bone conduction data can also be entered. These features can be turned on or off in Setup > Test so as to show only the required information.



Enter or change Test Data

**ASPARTATE Garry** Save Cancel

Date of Test \* 14/12/2017 DOB: 06/01/1955 Age 62  
 Next Test due 18/04/2011 Years of employment 9  
 Next Training 18/04/2013 Type of Test Routine test

\* = compulsory fields  
 This is the current Reference Test ☐  
 Visual examination of the left ear Normal  
 Visual examination of the right ear Normal  
 Test # 1  
 Occupation ENGINEER  
 Department INFORMATION SYSTEMS

**Audiometer & Tester**  
 Audiometer Audiometer 1 \* Serial No 1234 Calibration date 3/04/2000  
 Tester Goody Gumdrop \* Approval 1234

**Vehicle driver test**  
 Left 22.5 dB  
 Right 12.5 dB  
 Average at .5, 1, 2, & 3K  
 % Loss age corrected, NAL Tables 1998  
 Recalculate

**Percentage loss**  
 Left 6.8 %  
 Right 0.0 %  
 Bin 1.6 %

**Hearing protection**  
 Protection worn before test ☒  
 Protection usually worn ☒  
 Training provided ☒  
**Protector currently fitted**  
 Protector Type QB2  
 Protector Class 2  
 Select Protector  
**Protector required**  
 Protector Class 0

**OSC930** Help ?  
 Capture Test Capture Re-test

**Air Conduction thresholds**

	.5k	1k	1.5k	2k	3k	4k	6k	8k	Hz
Left	15	25	25	25	25	45	55	50	dB
Right	10	15	15	15	10	15	20	30	dB
Retest					20	45	50		
Left average	40.0				Mean	22	45	52	Left
Right average	13.3				Mean	10	12	17	Right

Left test is abnormal Right test is normal

**Hearing Test** **Spirometry**

Mouse draw controls  
 Left Right Copy  
 Air L R C  
 Off

**Medical Information**  
☐ Ringing ears  
☐ Frequent Colds  
☐ Discharges  
☐ Noisy Hobbies  
☐ Power Tools  
☐ Dizziness  
☐ Referred medically  
☐ Any change in medical history  
☐ Currently has ear ache?  
☐ Exposure to explosive/gun fire

Period of quiet before test

Select Comment Comments ABC Spell check

Pressing on the Add or Edit button on the Test tab opens the Test Add / Edit window

Note that the Australian Standard has been enabled, so you see second line of thresholds at 3, 4 and 6kHz and the resulting results are the yellow thresholds.

The Australian Standard can be enabled in **Setup > Australian Standard**.

The fields in yellow are compulsory, they are required to qualify the test.

In this view the Australian Standard Re-test values are shown at 3, 4 and 6kHz below the original test.

The values in yellow become the test results. This is only one of the options, you can also display the standard view which does not include the Re-test values.

The “**Capture Test**” button allow data from the audiometer to be captured automatically.

The audiogram can be drawn manually using the mouse or keyed in manually.

In the Data Entry window, click on **Print, Letter**. A window will open for you to select an existing file. The file must already exist. You can use an existing file as the basis for a new letter and alter the text.

Double click on the file name or click on the file name and click on **Open**. Your original file will be available with all the Mail Merge fields replaced. You can alter the text if you want before printing the letter. If you have set a **number of copies** in the **Setup, Page Format**, several letters will be printed.



### **How do I Mail the Letters**

The letters are designed to be mailed in a window envelope, this saves having to print an address label.

The top and left margins can be adjusted so that the address shows in the window of the envelope.

<p><b>Standard window envelope format</b> : Envelope size = 230mmX100mm, window is 58mm from top and 32mm from left side.</p>
---



## CHAPTER 4 - Starting Testing Employees

### 4.1 Setting up Defaults to Simplify your Work

---

If you know what company you will be working with that day, you should set that company as your default.

If the company exists already in your database, visit **File > Set User Default** (from the main menu) Then select the company, the tested and the audiometer you will be using that day as defaults.

If the company does not yet exists, go to **Setup > Company** and enter a new company. You can set it up as a default right there.

You can do the same for the Department, Tester and audiometer.  
Setting these up, will not only simplify your work but increase accuracy.

### 4.2 What Happens to the Defaults

---

When you enter a new employee, the current defaults are automatically filled into the relevant fields for you. These fields are  
Company and Department

When you enter a new test, the department and occupation are automatically copied to the test and the audiometer, date calibrated and tester name are automatically entered.

### 4.3 A good habit

---

Before staring work, **Go to File > Set user Default** and select your defaults for the day.

Remember when entering the company that a company is defined by the name AND the address. So make sure to enter the address as well as the name.

For example

Company: Qantas      Address: Sydney

The program will treat this company as different from

Company : Qantas      Address: Melbourne



## CHAPTER 5 - Things to Know

Windows programs are designed to be very graphic. The mouse becomes an extension of your hand to perform actions on the screen that your fingers would have performed in reality. Pushing a button becomes clicking with the mouse, on the representation of a button on the screen.

In this way, Windows programs are easier to learn and usually more intuitive than DOS programs were, but they require the use of both the keyboard and the mouse.

Each time you start the Workplace Audiometrics, the first screen shows the program's name and the credits as well as the working directory. A number of buttons appear on the top bar, they are called “Speed buttons”, they are shortcuts to the main functions of the program. These functions are also available from the **File** menu.



Entering or changing data is always done in the Data entry screen. Looking at data is done in the windows controlled by any of the 3 buttons. Note that some data can be changed from the browse windows, but nothing can be added.

### 5.1 Selecting From a Combo Box


---



Click on the button on the right (the one with an arrow down). A window will open, with a number of choices. Click with the mouse to highlight the desired choice, then click on the choice. The content of this line will become the content of the Combo Box. If you wish you can also type directly into the Edit window (the area where “OHS Employee” appears above).

### 5.2 Selecting from a Tabbed Book

---




In the above image the currently selected tab page is the Employee details, it is highlighted in Yellow on Black. To move to the tab page named **Medical History**, click on the Tab name “Medical History” with the mouse and that page will come to the foreground. Tabbed Books are used to display more information in a smaller area.

### 5.3 Selecting from a Date Field



To enter a date you can just type the date in the window making sure you enter the year as 4 digits (Example 2000 ). You can also click on the button (see above), when you do, a window with a calendar will pop up. You **must** click on the day to change the date when the calendar is showing. Click on >> to increase years, on > to increase month. Then click on the day to select it.

### 5.4 Moving and Selecting in a Grid




The current record is highlighted in red in the grid. To move from one record to another in a grid, you can click on the arrow on the top right to show hidden records above the top one, or on the arrow at the bottom right to show hidden records below the bottom. If the record you want is on the screen, you can just click on it to make it the current record.

To select from a selection window, you can double click on the red line, or you can highlight a line by clicking on it with the mouse, then click on the “**Select**” button.


### 5.5 Adding an Employee

Click on the **Add** button  
To change (or Edit) an Employee’s data,  
click on Edit



### 5.6 Changing from Employee to Test

Click on the **Test** or Employee tab to  
change from one area to the other.



### 5.7 Adding a test to an Employee

You can **Add** or **Edit** (change) an Employee’s test by clicking on the relevant button after selecting the **Test** tab.



Browsing employees. On the right, you can the tests for the employee highlighted in red on the left.

The Browse Test screen. The view can be set to show only tests performed within a date range and Normal / Abnormal tests.




### The Import Mode

In this mode, the software will automatically add each record to your data, in alphabetical order. If an Employee already exists (based on same Surname, Name and Date of Birth) only the test will be added. If the Employee does not exist in your data (or has a different ID number), both the Employee and the Test are added.



### Update Mode

In the update mode, only the fields in the update file containing data are used for




updating your data. If a field is empty in the update file, and the corresponding field in your data contains data, that data will not be altered.

**Important:** Be extremely careful with the files you import into your database. It may be hard to chase up and remove wrongly imported data.  
If the imported file is not in the right format, garbage will be entered into your database and most likely damage it.

**Important:** It is recommended to do a backup BEFORE the import or update.

## CHAPTER 6 - The File Menu

Set User Defaults		
Data Entry		Ctrl+D
Browse Employees		Ctrl+E
Browse Tests		
Import from a WA Text File (ASCII)		
 Backup		Ctrl+B
Data File Statistics		
Exit		

### 6.1 Data Entry

This menu option is replicated as a Speed Button. This is where all the data entry is performed.

### 6.2 Browse Employees

This menu option is replicated as a Speed Button. It is a way to browse through your employees. This option is not meant to be used for data entry, but data can be modified there. Employees are shown in Alphabetical or Work number order at your choice.

### 6.3 Browse Tests

This menu option is replicated as a Speed Button. It is a way to browse through your employee's Test. This option is not meant to be used for data entry, but data can be modified there. Tests are shown in entry date order. The display can be set to only show All, Normal or Abnormal test and by a date range.

To change the range, enter a date from and a date to, then click on **Search**. To view all tests click on **Reset**. This window show all tests by default.

### 6.4 Set User Defaults

Allows you to setup the default company, the default tester and the default audiometer. This information is copied automatically to the Employee and test records, so it is worth while making sure they are correct at the start of the testing day.

### 6.5 Import from a File WA Text File

Data created by another computer in the correct format and stored in a file can be imported (Added) to your list or used to update it. This File could have been created by your Mainframe, another PC or by the **List, Export Data to a WA Text File (ASCII)** menu option in another copy of this software. See the end of this manual for layout of data to be imported in case you want to import an ASCII delimited file from another source. Data can be imported in two modes.

### 6.6 Run Another Program

This feature is designed to allow running of another software, or a software module from within this software. To activate this function, go to **Setup, Preference**.

Under **Run another program 1 (or 2)**, enter the name of the software and the title to appear in the **Files** menu.

for example:                      Command line: C:\DUMMY\TEST.EXE  
   Title:                      Testing Software

Of course this software must be present in the correct directory before this function can operate properly.

## 6.7 Backup Data

---

The data in your database is very precious and as such should be protected against loss. The hard disk is the most fragile component in your computer, yet, it stores all your data.

The data should be backed up from the hard disk regularly. This is done in a compressed format called ZIP. The easiest way to backup is to a USB Stick. Best way is to get 5 USB Sticks, label then Monday, Tuesday... and backup on the one for that day.

It is **STRONGLY recommended** to use **multiple sets** of backup for added safety. You can even store a spare set in another location to protect against fire and theft.



### When to Backup

The backup should be done every day that data is added, on a different USB memory stick.

:



The backup should be done every day that data is added, on an alternate USB Stick. You can do the backup even while the program is used by an other user.

## 6.8 Doing the Backup

---

To do the backup, insert a disk in the drive, then select the backup option. The screen will list the files as they are backed up. Note that the letters will also be backed up.



### Where is the Backup saved?

You can select where the data is saved by default in **Setup > Preference**, in the **Option** box. It normally is the A: drive and can be changed to any drive.

To backup to a USB Stick, just insert it in your computer USB connector, watch to see a new window



open, look at the top of the window to find out what drive the USB Stick is, then go to **Setup > Preference**, in the **Option** box select that drive as the backup drive. Backup will then automatically be sent to the stick for you.



### **The Advanced Option**

In some cases, for example if you have a network or a zip drive, you may want the backup files to end up in the sub-directory of a drive, not in the root directory. This can be particularly convenient when doing a backup every day in a different sub-directory. for example : F:\MONDAY, F:\TUESDAY etc.



### **What is the format of the Data**

The data is in standard ZIP format, so it can be un-compressed by WinZip for example. But WinZip will not know where all the files should go, so it is safer to use the **Restore** function, in the **Utility** menu to handle a crisis unless you cannot re-start the program to do the Restore. See the Restore function later in the manual for more information on this subject.



### **Can I restore from a Backup from an Older Version of the program**

No, you cannot because the data is in a different format and data files have different names purposely.

## **6.9 Data File Statistics**

---

Shows the number of entries in each of your data files. The numbers include the deleted records. The deleted records will not show in the reports.



## **CHAPTER 7 - Data Entry Reports**

### **7.1 The Employee's Audiogram**

---

The current employee, in the data entry screen, can be printed by pressing the Print Button and then the Audiogram button. Before pressing the Audiogram button, you can select to send the audiogram to a preview screen, to the Windows default printer or to save the Audiogram as a pdf document. The print to pdf will create a file namer made of the employee's name.

If you pressed the Print button while looking at the employee tab, the last test will be the one automatically included in the audiogram

If you pressed the Print Button while looking at one of the employee's tests, that test will be included into the audiogram. This mean that you can re-print any test, by simply highlighting it in the left grid.

### **7.2 The Form 411 Report (WA)**

---

This report is used by Workcover in WA. It operates the same ways as the above audiogram, it just look different.

### **7.3 The CX Report**

---

Similar information as the Form 411 report but for general use, not just Western Australia.

### **7.4 The Report Card**

---

This report combines all the tests for the employee currently selected. At the end of the report is a comparison between the oldest test (baseline) and the most recent test, showing the loss in threshold between the two.

Like the other reports, can be printed, previewed or saved as a pdf document.

### **7.5 The Spirometry Report**

---

This report is only available if you have ticked "Spirometry" in Setup > Test > Spiro, Tymp, Bone > Spirometry.

Just like the Report card it shows all the spirometry tests for the employee currently selected.

### **7.6 The Percentage Loss Report**

---

This report will print a detailed Percentage loss report with a breakdown of losses per frequency. It is available in Standard (% Loss Std) and Insurance format (% Loss Ins). This report can use either the Air conduction or Bone conduction thresholds. The default is Air conduction.





A quick evaluation of a particular person, based on the Hearing test Baseline and the current test can be obtained at a touch of a button in the Data entry screen.

**Hearing Test Assessment -**

### Hearing test assessment report for ASPARTATE Garry

**Australian Standards AS1269.4:2014**

- (a) Shift from Baseline averaged at 3000, 4000 and 6000Hz  $\geq$  5 dB
- (b) Shift from Baseline averaged at 3000 and 4000  $\geq$  10 dB
- (c) Shift from Baseline averaged at 6000  $\geq$  15 dB
- (d) Shift from Baseline at 500, 1000, 1500 or 2000  $\geq$  15 dB
- (e) Shift from Baseline at 8000  $\geq$  20 dB

**Victorian regulation 2017 3.2 & New Zealand Standards**

- (v) Shift  $\geq$  15db at 3000, 4000 or 6000Hz between last test and tests within a 2 year period

**Rio Tinto NIHL**

- (f) Current thresholds averaged over 1000, 2000 and 3000 Hz  $\Rightarrow$  25 dB in both ears, age corrected

**OSHA STS (01/01/2003) age corrected**

- (1) Shift from Baseline averaged at 2000, 3000 and 4000Hz  $\geq$  10 dB, age corrected
- (2) Average absolute at 2000, 3000 and 4000  $\Rightarrow$  25 dB same ear as in (1)

**OSHA STS (01/01/2003) non-age corrected**

- (h) as above but not age corrected

**OSHA STS (Canada)**

- (3) Shift from Baseline at 2000, 3000 and 4000Hz  $\geq$  30 dB, age adjusted

**Commercial vehicle and Train Drivers**

Pass = last test average at 500, 1000, 2000 and 3000Hz in better ear < 40 dB

**Tram Drivers**

Pass = last test average at 500, 1000 and 2000Hz in better ear < 40 dB

**Results**

Australian Standard					Reg. 2017 (v)	Rio OSHA		OSHA S.T.S			Age corrected S.T.S		NIHL		Phone Call centre Operators	Drivers Commercial vehicles Trains	School Crossing guard
(a)	(b)	(c)	(d)	(e)	(v)	(f)	(h)	(1)	(2)	(3)	Left	Right	Left	Right			
X											-3	5	175	75	Fail	Pass	Pass

**NAL 80 (Waugh & Macrae)**

(1)	(2)	(3)	(5)

L=Left, R=Right, B=Both, Y=Yes

**Information:** The software has assessed this employee according to the standards outlined on the left. The result of the assessment is shown below. An X in a box means that a person has a hearing shift exceeding the level permitted by that standard, therefore testing positive. Note: Only the Australian standard and the OSHA STS are using the Baseline test. The Baseline test is defined as the latest test with a test type = B, or if not present, the oldest test.

Last tested on: 14/12/2017  
Baseline test date: 01/02/2008  
Number of tests: 2

[Read help for complete explanation of the standards](#)

[Help ?](#)

[Save](#) [Print](#) [Close](#)

☐ CSV ☒ RTF

## Australian Standard Report

(a) a shift in mean threshold at 3, 4 and 6 KHz  $\geq$  5 dB (b) a shift in mean threshold at 3 and 4 KHz  $\geq$  10 dB  
(c) a shift in thresholds at 6 KHz  $\geq$  15 dB (d) a shift in mean threshold at 0.5, 1, 1.5 and 2 KHz  $\geq$  15dB  
(e) a shift in thresholds at 8 KHz  $\geq$  20 dB The shift is computed from Reference test to latest test

Selection: All Employees

Name	Work No	Sex	Age	(a)	(b)	(c)	(d)	(e)	Baseline	Last test	No tests
ASPARTATE Garry		M	54	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	01/02/2008	18/04/2009	2
BELLA DONNA		F	89	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	11/11/2013	11/11/2013	3
BLOB Fred		M	53	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	01/08/2005	01/08/2006	2
STRADIVARIUS Viola	12345	F	38	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	18/02/2002	10/10/2013	6

Number listed: 4

Total(a): 1

Total(c): 1


Total(e): 2

Printed: 22/01/2014

Total(b): 1

Total(d): 1

The Australian Standard Report conveniently shows the positive test for a group of employees.



The report menu allow printing personal Employee reports in batches, for example at the end of a testing day as well as report on groups, or particular report including employees meeting a set of particular criteria, like company, department, test abnormal etc..

An extensive set of statistics report is also available.

Just in case none of the reports are suitable for your purpose, you can create you own report based on your criteria, and save your criteria for future use. Multiple criteria can be saved that way.

The potential liability report will expose the financial risk from employees claims in the unrealistic case where they would all claim at the same time. This is used mostly by self insured company to know how much to put aside to cover themselves.

## CHAPTER 8 - The Report Menu

### 8.1 Batch Audiograms

---

After accepting your selection, this function will print the latest Audiogram of the Employees matching your criteria except if the selection includes the “BY DATE” option. In this case all Audiograms matching the date are included regardless of the Employee.

The Batch Audiogram function ignores the deleted records.

To print Audiograms of a previous Test for one specific Employee, go to the **Data entry**, find the Employee, click on the **Test** tab and press **Print, Audiogram**.

### 8.2 Batch Report Cards

---

This function will scan through the database and find the Employees matching your selection.

For each Employee it will collect all the audiometric data and print it in the form of a report card. There are 2 types of Report Cards, one for Audiometry and one for Spirometry.

### 8.3 Batch Form 411 (WA)

---

This function will scan through the database and find the Employees matching your selection and print a Form 27 (WA) for each of them.

### 8.4 Batch CX Report

---

Same as above but for general use.

### 8.5 Victorian Regulation 2007 3.2.11 (Test shift)

---

This function will create a complete report, ready to send to the appropriate government department (in Victoria) meeting the requirements of Regulation 1992:

- Details for the current test are completed if the results at 4000Hz equal or exceed 25dB at 30 years or less, 35dB at 45 years or less, or 50 dB at any age.

This report highlights the Employee’s deterioration in hearing from one year to another, and can be used regardless of the regulations in your state.

### 8.6 Australian Standard Report

---

This function will create a report with check boxes. Each check box represents one condition requiring that retest be performed on that employee. The conditions include the Rio Tinto (USA) parameter. The conditions are as follow:

- a) Shift from baseline averages at 3000, 4000 and 6000Hz  $\geq$  5dB
- b) Shift from baseline averages at 3000 and 4000  $\geq$  10dB
- c) Shift from baseline averages at 6000Hz  $\geq$  15dB

- d) Shift from baseline averages at 500, 1000, 1500 or 3000Hz  $\geq 15$ dB
- e) Shift from baseline averages at 8000Hz  $\geq 20$ dB
- f) Current thresholds averaged over 1000, 2000 and 3000Hz  $\geq 25$ dB
- 1) Shift from baseline averages at 2000, 3000 and 4000Hz  $\geq 10$ dB age corrected
- 2) Average absolute at 2000, 3000 and 4000Hz  $\geq 25$ dB same ear as in 1

In addition, a button called “Assessment” is available in the Data entry window to provide the same information for the currently selected employee.

## 8.7 The Assessment Screen

The option is located in the Data entry window. It provides a quick check of the currently selected employee.

Additionally to the Australian Standard test, 5 more tests are performed:

- Call centre operators
- Commercial vehicle drivers and Train drivers
- OSHA regulation for USA based companies.

These results do not appear in the Australian standards report. They are part of a dedicated OSHA report available from the Report menu. Both report can be exported in pdf or text format. The data conforms to the CSV standard and the data lay out is as follows:

Surname, Name, Work number, Code, Sex, (a), (b), ©, (d), (e), (f), (16,1), (1), (2), Occupation, Department, Shift, Date of last

Australian Standard					Reg. 2017	Rio OSHA		OSHA S.T.S			Age corrected S.T.S		NIHL		Phone	Drivers		School
(a)	(b)	(c)	(d)	(e)	(v)	(f)	(h)	(1)	(2)	(3)	Left	Right	Left	Right	Call centre Operators	Commercial vehicles	Trams	Crossing guard
X											-3	5	175	75	Fail	Pass	Pass	Pass

## 8.8 Base Line Shift Report

This function will create a report with 2 lines per entry. The first line is the Current test (most recent test stored in the list), the second line is the Base Line test (oldest test stored into the list).

An Employee will only appear in the report if their Current test is abnormal. The second line is completed if there is a difference greater than 15Db at any frequency from 500Hz to 6000Hz, between the Current test and the Base Line test.

The Baseline shift is an excellent report to show the highlight the Employee’s deterioration in hearing from the first test (baseline) to the last test.

## 8.9 Whole Person Impairment Report (Workcover NSW)

This report is part of the Australian NSW regulation and complements the Potential Liability

report.

### **8.10 Whole Person Impairment Report (Comcare)**

---

This report is part of the Australian Commonwealth regulation and complements the Potential Liability report.

### **8.11 OSHA STS Report**

---

This report is used in the US. Some companies with a head office in the US may want to obtain this report from your organisation for comparison purposes. A separate version is available for Canada.

### **8.12 Category Report**

---

This report breaks the employee list into 4 category based on the employees hearing thresholds averaged at 500, 1000, 2000 and 4000Hz.

### **8.13 Design Your Own Reporting Criteria**

---

This menu option lets you design your own reporting criteria. Once created you can save the criteria to be able to use it again. You can save any number of criteria. This option is very flexible.

### **8.14 Assessment Report**

---

This report will print an assessment of the employee selection based on the standards available.

### **8.15 Commercial Vehicle and Train Driver Reports**

---

This report will give a pass / fail report based on the selected criteria.

### **8.16 Report Creation**

---

Most reports can be very lengthy documents. It may be advantageous to know in advance the number of Employees included and the number of report pages required. To do this simply print the report to the screen, click on the last page icon at the top to see the number of Employees included and the number of pages.

### **8.17 Selecting Data for Reports**

---

The selection process scans the Employees and their **LAST** test to find the tests falling within the date range (month or dates). The data is then analysed to assess if it should be included in the report and in what Regulation.

### **8.18 Employee Status Statistics**

---

Will calculate the number of Normals, Abnormals and Totals for each Status classification and show the results with figures and bar graphs.

### **8.19 Age Statistics**

---

Will give a breakdown of the Normal and Abnormal Employees based on age in the ranges 0-20, 20-30, 30-40, 40-50, 50-60 and above 60.

---

## 8.20 Employee Medical Statistics

---

Will give a compilation of the answers to the Medical History questions in the Employee record. This would reflect the Employee's condition when joining.

---

## 8.21 Test Parameters & Medical Statistics

---

Will give a compilation of the answers to the Medical History questions in the Employee's test record. This would reflect the Employee's condition at the time of the last test.

---

## 8.22 Potential Liability

---

The potential liability report is designed to give an estimate of the amount of compensation that would have to be paid in compensation, if all your employees made a successful claim at the same time. This figure is, of course, unrealistic but it is a good approximation of maximum risk. It could be used for financial planning, especially in case your company is self insured.

The report can be run for a specific population of employees. It can be run in two different ways:

- 1- using the relevant figures for your state for all your employees, regardless of which state they are from.

- 2- using the relevant figures for one state and scanning only the Employees from that state.

There are two formats: Global figure and a detail list with each Employee.

Because of the confidentiality of this report it is protected by a Master password (if you have activated the password function).



### **Single Formula for All your Employees**

In this case it does not matter if your employee's data includes the state they live in.

- Select the state your company is located in
- Make a selection
- Run the report

The report will scan ALL your employees to compile the report.



### **Specific Formula for each State**

Each employee's record must contain the state they live in.

- Select one state where some of your Employees live
- Make a selection
- Run the report

The report will only scan the Employees from that state. Select the next state and run the report again etc...

If you want to include New Zealand, since there is no State, always have "Scan all your employees" selected, or enter NZ as the state.



**List Employee data**

### List Employees

**Option 1**

- ☒ All
- ☐ By date of last test
- ☒ By Shift
- ☒ By Company
- ☒ By Department
- ☒ By Company & Department
- ☐ By Occupation
- ☐ By Employee's Postcode
- ☒ By tag
- ☒ By Workg

• = High speed search

**Option 2**

- ☒ All
- ☐ By date of last test
- ☐ If last Test is Abnormal
- ☐ If last Test is Normal
- ☐ If % loss both sides greater than X
- ☐ By a keyword in the Test comments
- ☒ By Employee's Status
- ☒ By Shift
- ☒ By Workg
- ☒ By Cost Centre
- ☐ Claim submitted
- ☐ Claim settled

**Output**

- ☒ Screen
- ☐ Printer
- ☐ pdf File
- ☐ Labels
- ☐ Letters
- ☐ Mail Merge

**Format**

- ☒ Audio data
- ☐ Personal data
- ☐ Company data
- ☐ Short personal

**Include**

- ☐ Deceased
- ☐ Terminated

Printer Setup Page Setup Top margin = 10 Left margin = 10 Orientation = Landscape **Start** **Exit**

To create a report you need to select the criteria of your choice from the box on the left then the box on the right.

You can then select an output (for example to the printer)

All formats show the same list of Employees based on you choice, but they show a different portion of the Employees information.







## CHAPTER 9 - The List Menu

### 9.1 List Employees

---

Displays or prints a List or subset of your Employee List using your selection criteria. Four formats are available to cover most requirements. For additional requirements you can export the data to Excel and pick your choice of fields.

**Warning:** Some lists may take a while to print and the computer may appear to have locked up. This is most likely not the case, you will just have to wait until data is found.

In many cases a progress bar is visible, but it only appears after your query has been implemented and is shown during the composing of the actual report.



#### **Page Setup**

This button accesses the same menu option as **Setup, Page Setup** on the main window. It allows the customisation of all the functions found in the List Menu. Spend some time to explore each tab to see what can be done.

**Warning:** Pay special attention to the size of the margins. If they are too wide and push the data out of the page, the page will be printed blank.



#### **Printer Setup**

The lists and reports can only print to the Default printer as defined in the **Control Panel** of Window. This is a limitation which cannot be overcome at this point. So even if you select another printer, data will still go to the Default printer.

### 9.2 Listing Output Format

---

Listings can be sent to the Screen as a Preview, then printed if you wish, or sent directly to the printer.

Three formats are available:

- 1- Audio Data
- 2- Personal Data
- 3- Company Data

Additionally the listings can have grid lines or not, and include Deceased and Terminated Employees, or not.

### 9.3 Labels

---

This software is designed to use A4 sheets of sticky labels as available from Unistat, Avery and some other manufacturers.

Labels are defined in the **Page Setup**. Four pre-defined formats are available. You can also define your own by selecting the Custom label and entering the relevant figures in the label data window.

Labels print from top to bottom and from left to right.



You should now have a blank page, click on the **Insert** menu and select **Field**

- In the left window click on **MailMerge**, in the right window, click on **MergeField**
- The word MERGEFIELD appears in the **Description** window, click in the window after the word MERGEFIELD and type **NAME** , click on OK.

You should now have the word <<NAME>> in the document. At the time of printing this will be replaced by the Employee name.

Following the same procedure define each field with the same names you have selected in the software in **Setup, Page Setup, Mail Merge** .  
Locate them where you want them and add the body of text, of the letter. Save the file as MYLETTER.DOC for example.

If you forgot the name of your fields, just save to document, (or Alt Tab to the software) go back to the software and check the name of the fields in Setup, Page Setup, Mailmerge, then come back to the document.



### **Printing The Mail Merge Letter**

With the main document we worked on earlier showing on the screen, click on the **Tools** menu, select **MailMerge** .

We are now at stage 2, click on the **Get Data** button.

- Click on **Open data source**
- In the **Look in** window click on C, then double click on AUDIO
- At the bottom, in **File type**, click down to **Text file (.TXT)**
- To choose the file you just created with the software (MYDATA>TXT for example ) double click on the filename

We are now at Stage 3, click on the **Merge** button next to Stage 3

- Press on the **Merge** button at the top right

The creation process will take a while.....

When it is completed, you will have a file called **Form letter1.doc** containing all your letters.

- Click on **File > Print** and print the letters.

## **9.6 Printing Mail Merge Envelopes with MS Word**

---

You would be better off creating your letters with the address in the right place for a window envelope as envelopes do not often feed well in the printer, but for your information here is how to:

This example assumes you have already printed the letters and the Data source has been created as mentioned above.

- Open Microsoft Word, select **File > Open** - Open C:\AUDIO\MYLETTER.DOC
- Select **Tools > MailMerge**, press the **Create** button and select **Envelope**
- Press the **New Main Document** button
- Press the **Get data** button and select **Open data source**
- In the **Look in** window click on C: drive, then double click on AUDIO
- At the bottom, in **File type**, click down to **Text file (.TXT)**
- Choose the file you just created with the software (MYDATA.TXT for example ) and double click on the filename
- Press the **Setup Main document** button
- Press on the **Merge** button at the bottom

- Select the envelope type / size then press the **OK** button
  - Press the **Insert Merge fields** button and select each fields to create the address label for the envelope, making sure to leave a space between Name and Surname
  - Click the **OK** button when completed
  - Click on the **Merge** button at the bottom
  - Click on the **Merge** button at the top right
- The creation process will take a while.....
- When it is completed, you will have a file called **Envelopes1** containing all your envelopes.
- Click on **File > Print** and print the envelopes

Tip: if you use window envelopes and set the address in the right place you can remove the need to print the address labels on the envelopes.

---

## 9.7 List Tests

Displays or prints a subset of your Employees test list using your selection criteria.

This option is ideal to list the all the tests done in a date range, or for a Department or a company.

The Tests are not stored in the same data file as the Employees. They are linked by the software. It is possible that Test belonging to a deleted Employee are still present. In this case the Employee name will not appear on the report.

You can use the **Utility > Data Integrity** menu option to find and delete the tests belonging to deleted Employees.

---

## 9.8 Export Data to a WA Text File

This function will send your data to a file in a Text only format called ASCII delimited. Each field is enclosed in double quotes and separated by a comma. This file may be used by other software to import the data, or it can be used to import into another copy of the Workplace Audiometrics software.

The content of this file can be controlled by a selection so that only the Employees you want transferred are exported.

This function can be used to transfer between two locations both running the Workplace Audiometrics software.

A choice can be made to export only the Employees and their last test, or the employees and each of their tests. This will inflate the file considerably since the entire employee's record has to be present with each of the Employee's Tests.

When importing, the program checks for the existence of the Surname, Name and Date of Birth. If a match is found for the Employee, only the test data is added (if a test date is present). If the date of test is missing, only the employee will be added to your Database.

---

## 9.9 Export to Excel

The Employee and Test can be exported to Excel in **List > Export to Excel**. You will be

prompted to enter a file name and a location. The file will be saved as a csv file.

To open it in Excel, just select the file type .csv, you can then save it as an excel workbook if you wish. The same selection criteria as for other list apply.

## 9.10 Re-test Management

---

Each person should be tested once every period, according to their Employee status (exposed to noise or not). The period length for the Employee's status type can be setup in the **Setup, Test , Test Setup** tab, **Re-Test** box. This is normally a year or two for noise exposed Employees.

To facilitate the re-test process, a re-test or Next test date is stored as part of the test data. When the test is saved this date is copied to the Employee data automatically.

The software calculates the re-test date if the Employee **Actions Allowed** includes Re-Test.

If the test is a baseline test, a separate re-test date gap can be also entered, so that the software can allocate for example

- Date of test + 3 months if it is a baseline test
- Date of test + 2 years if it is a normal test

The Re-test Management uses the "Next Test Due on" date in the Employee to list the Employees due for re-test.

The list can be broken into departments and shift etc. so that each supervisor can be notified separately.



### **Labels for Test Reminder Mailings**

Labels can be printed by selecting the **Labels option**.



### **Letters for Test Reminder Mailings**

Use the MailMerge option to create the dataset needed by MS Word and do a Mailmerge following MailMerge instructions in previous chapter.

## 9.11 Re-test Scheduling by Departments in Company

---

This option will allow you to select one company and to list the employees dues for retest by departments with employees within a department in Alphabetical order.

## 9.12 List Companies to Contact For Retest

---

This menu option will list all the companies when one of their employee needs to be re-tested based on the date of their last test.

This feature allows you to know which company to check for re-test without having the list them all.

## 9.13 List Departments to Contact For Retest

---

This menu option will list all the departments within a company when one of their employee needs to be re-tested based on the date of their last test.

This feature allows you to know which department to check for re-test without having the list them all.

## 9.14 Re-Train Management

---

When Employees are tested, they should be trained in the use of Hearing Protectors. Some States (Vic) actually make it compulsory to record that the training has been done and to do it regularly.

The software allows you to manage this by recording that training has been provided and by setting a re-train date. The re-train gap can be set in **Setup, Employee, Status, Tags & Re-train** tab, **Re-train Employee** box.

## 9.15 Hearing Protection Management

---

In the Employee's workgroup there is a figure (Exposure Class) representing the Employee's exposure to noise. In the Employee's Test there is figure representing the amount of protection provided by their hearing protector (Hearing Protection Class).

This function will list the Employee when the two figures are not equal.

In the Employee tab, in the Data Entry window, the protector class appears in red if the two figures are not equal.



## CHAPTER 10 - The Utility Menu

### 10.1 Re-index File


It is possible that due to power surges, from time to time, an index can get out of order. The symptoms are: Employees or tests appear missing, but you know they should be there. In this case, just re-index the file and everything will return to normal.

### 10.2 Untag Employees

This function will remove specific tags that may have been attributed to an Employee. All the tags of the type you specify will be removed. During the process, the names of the tag holder's will be displayed on the screen

### 10.3 List Duplicates

You can list the duplicates and select the basis on which the program considers two records to be duplicates. The obvious choice is Surname and Name, but you can add the date of birth and the Employee number.



### 10.4 List Deleted Records

Records can be marked for deletion. They still appear in the Browsers (with a dark blue first field) but do not appear in any of the reports or lists. These deleted records are still in the list. They are just called deleted.

You can list them here before destroying them to make sure none was deleted by mistake.

### 10.5 Purge Deleted Records

To remove the records marked for deletion, go to the **Utility > Purge deleted**. This will remove them physically from the database. If you want to can purge selectively each list by using the check boxes.

### 10.6 Rename a Company / Department / Occupation

This function will:

- change a Company / department / Occupation or Workgroup name from the "Department" file from old to new
- look for all the Employees who belong to that department
- when found, it will change the name in their personal record.

- display the name of the Employees updated
- print the total changes

---

## 10.7 Update Department & Company Numbers

---

When you select a Company or a Departments, a number is incremented to reflect the total number selected. When you rename a Company or Department or delete Employees, this number is no longer accurate. This option will re-calculate the number and delete the Company or Department with number of Employees = zero.

---

## 10.8 Update Employee Retest Date

---

This function is available in case the update from an earlier version of the program did not transfer the Re-test date from the Test to the Employee's record.

---

## 10.9 Recalculate Percentage Loss

---

This function will recalculate the percentage loss for a selected group of employees. The Formula used can be selected from 1976, 1988 or 1988 extended calculation formula.

---

## 10.10 Switch a test from One Employee to Another

---

Sometimes a person may be duplicated in the list under a misspelt or slightly different name. This feature will allow you to move the test from the wrong name back to the right name. The wrong name can then be deleted.

---

## 10.11 Restore Data

---

Allows the data stored on a backup disk to be restored into the computer. Extreme care should be exercised here because loss of current data will occur irreversibly. This function will re-create the data as it was when it was backed up. Of course the current data is erased in the process.

Please think twice before using this function and make sure that the current data is also backed up, just in case.

The RESTORE function complements the BACKUP function. Do not use it lightly: it erases all of the current data and replaces it with the data stored in the backup.

It is advisable to contact a computer expert before doing a "Restore". The RESTORE function, like its counterpart the BACKUP function, uses the ZIP format used by WinZip.

If you are not able to start the program because the data files have been damaged, run WinZip and unzip the backup into your working directory (normally C:\AUDIO) then copy the file SCREENW.INI to your C:\WINDOWS directory to restore your local settings.





## CHAPTER 11 The Setup Menu

This menu provides a number of choices allowing you to customise many of the program's parameters to your own requirements.

### 11.1 Setting up a Password

---

The software can be protected by passwords.



#### **Cancelling the Password**

If you decide not to use the password protection, un-tick Enable Password Security, no passwords are required or requested.




#### **Turning Password On**

If you decide to setup the password protection, a new menu item will appear in the Setup menu, called Security. There you will be able to add users and password as well as users rights.

To enable Password Security, tick the checkbox.

You will immediately be taken to the Security Setup screen.



### 11.2 Adding Users to the Security System

---

Each user has 19 rights, allowing access to various parts and functionality of the program. Outdated users can be kept in the system and de-activated or simply deleted.

The Superuser cannot be altered. The admin, like the Superuser, has all rights enabled and only allows for the password to be changed. Make sure you keep that password safe for future access.

You can build any selection of access rights for specific users. When a menu option is not enable for that user, it will show as greyed out in the menu.

Only the Admin user is allowed to edit the tests past the day they are entered (to allow for error corrections).

To edit a test once the entry date is passed, you will need someone to enter the MASTER password for each test or to re-enter the program with the MASTER password.

Should you lose the MASTER password, you will have to call BIOTRONIC Pty. Ltd (02) 9398 8132 with proof of your authority as a MASTER user to get help logging in.

To change the user while in the program, go to **Files > Log on** and enter the new users log on name and password.

---

### 11.3 Entering Companies

Click on **Setup, Companies**. To add a new company, click on **Add**, to change an existing one, click on **Edit**. Your own or default should be entered as the **Default Company**. If you work with more than one company, then you can enter each one. In Victoria you will have to enter all of the company details, as they are used by the Regulation 1992 Summary Report. Click on **Exit**.

---

### 11.4 Entering Departments

Click on **Setup, Department**, there should be at least one Department. Departments can be added at the time they are selected, when entering an Employee's data. So there is no need to enter all of them in advance. The Number is updated automatically each time you select a Department. Click on **Exit when completed**.

---

### 11.5 Entering Occupations

As for departments and companies.

---

### 11.6 Entering Workgroups

A workgroup is an area of similar noise exposure in L Aeq, 8h. The noise level can be entered, so that when tests are entered, the wearing of an adequate hearing protector can be checked for that person. The software will automatically enter the corresponding noise class based on NAL and Australian Standard 1269, 1998. If the exposure is tonal (one specific frequency only) or above 105 dB this system is not applicable and a noise professional's help may be required. Click on **Exit**.

---

### 11.7 Hearing Protection Devices

The list supplied with the software was extracted from the National Acoustic Laboratory (NAL) booklet called "Attenuation and Use of Hearing Protectors" (eighth edition) and was included with NAL kind

**Browse Companies**

Search:

Doubleclick the Company you want to use

Company name	Tag	Company address	No	CONTACT	TEL
CUMMINS LAVERTON		BOUNDARY RD, LAVERTON, 3028FAS	5		
FISH AND SHIPS			4		
GANTAS			0		
TELSTRA		Norfolk STREET	52	Fred bloggs	9398 8132

Default Company: TELSTRA  
Norfolk STREET

Buttons: Set Default, Email, Add, Edit, Del / Undel, Print, Exit

## Setting up your list of Companies

Make sure to enter a default Company to speed up data entry.

If you have multiple Companies, enter the most commonly used one.

You can always setup the default before testing in File > Set User Defaults

**Browse Workgroups**

Search:

Workgroup Details

Workgroup Name	Tag	L Aeq, 8h	Tonal expos.	Class	#	Special
COOK		95		3	0	
PROCESS LINE		89		1	4	
TWRTR		0		0	1	

Default Workgroup: PROCESS LINE

Buttons: Add, Edit, Del / Undel, Print, Exit

**Workgroup Details**

The Classification concept and Hearing protector data have been kindly supplied by NAL.

You can enter a noise exposure value L Aeq, 8h, in dB(A). This value, if below 110 dB and not a tonal noise, will give this Workgroup a Protector classification that will be used to assess the protector required for this area.

This system is compliant with Australian standard 1269 and will allow you to assess the suitability of an employee's hearing protection by running the HFD report.

The HFD report will compare the protector class in the employee's workgroup with the class of the protector worn at the time of the last test

Beware that this system is only designed to cover non tonal noise levels below 110 dB(A). For higher noise levels or tonal noises, please refer to a noise professional.

## Setting up your Workgroups

Workgroups can be used in different ways. Here it is used to mean an area of same noise exposure level. This level can be compared with the Hearing Protectors worn by the Employees.

You can also use it for your own purpose.

**Browse Hearing Protectors**

Search for a Class:

Class: ☐ Make: ☐ Model: ☐

Doubleclick the Hearing Protector you want to use

Model	Tag	Class	Make	Type	Description	SLC 80
WORKCOM PRO			BILSOM	Comm headset	Communications headset	0
EMLF-48	1	PROTECTOR		Ear muffs	Safety Products EMLF-48	11
YELLOW 2500H MSA	1	UVEX		Helmet muffs	Yellow 2500H on MSA Hard hat	11
PHONSTOP	1	AM		Ear plugs	Phonstop	13
PROPP-O-PLAST	1	BILSOM		Ear plugs	Propp-O-Plast	12
CANALCAPS	1	CIGWELD		Ear plugs	Canalcaps	10
EAR DEFENDERS	1	MINE SAFETY		Ear plugs	Mine Safety Appliances Ear Defenders	12
DBA	1	RACAL SAFETY		Ear plugs	Racal Safety DBA	12
H.I.C.	1	TASCO		Ear plugs	H.I.C.	13
2470 COMPACT	2	BILSOM		Ear muffs	2470 Compact	17
SILENCER	2	CIGWELD		Ear muffs	Silencer 455200-1990	16
SILENTA ERGO II	2	KEMIRA		Ear muffs	Silenta Ergo II Cap attachment P/N 455308	15
SILENTA UNIC	2	KEMIRA		Ear muffs	Silenta Unicap Cap attachment P/N 455273	16
HUSHMUFF MK II	2	MINE SAFETY		Ear muffs	Mine Safety Appliances Hushmuff MK II 222775-09	17
NOISEFOE IV	2	MINE SAFETY		Ear muffs	Mine Safety Appliances Noisefoe MK IV	16
EML10 (PRE 1987)	2	PROTECTOR		Ear muffs	Safety Products EML10(pre-1987)	15
EMLU-47	2	PROTECTOR		Ear muffs	Safety Products EMLU-47	16
HUSHMUFF22215006	2	MINE SAFETY		Helmet muffs	Hushmuff 22215006	15
EMCC-50	2	PROTECTOR		Helmet muffs	Safety Products EMCC-50	14
EMM71	2	PROTECTOR		Helmet muffs	Safety Products EMM71	16
002.0881P	2	RACAL SAFETY		Helmet muffs	Racal Safety 002.0881P	17
YELLOW 2500H CIG	2	UVEX		Helmet muffs	Yellow 2500H on Cigweld Hat	15

Buttons: Add, Edit, Del / Undel, Print, Exit


The software will allow you to compare the Class of noise employees are exposed to with the class of the hearing protection device they are wearing.  
The Class of noise the employee is exposed to is stored with the Workgroup/Workgroup  
The Hearing protection device comparison mechanism in this software is based on The Australian Standard 1269.  
Warning: This system is only covers Non Tonal noise levels up to 110 dB. For higher noise levels or tonal noises, please refer to a noise professional.

## Setting up your hearing Protectors

The software comes with an existing list which can be added to.

The use of this option is not compulsory.

It is fine if you ignore it as it does not apply to every situation.



**Potential liability setup**

Select the State & formula parameters

	Minimum compensable	Maximum amount payable	Percentage of maximum amount payable
<b>Australia</b>			
<input type="radio"/> VIC	10.00 %	\$ 104990	65.00 %
<input checked="" type="radio"/> NSW	6.00 %	\$ 100000	65.00 %
<input type="radio"/> SA	5.00 %	\$ 105000	75.00 %
<input type="radio"/> WA	75.00 %	\$ 119000	119.00 %
<input type="radio"/> QLD	11.00 %	\$ 47995	100.00 %
<input type="radio"/> TAS	5.00 %	\$ 410.89 x 113	100.00 %
<input type="radio"/> NT	0.00 %	\$	0.00 %
<input type="radio"/> ACT	0.00 %	\$	0.00 %
<b>New Zealand</b>			
<input type="radio"/> NZ	0.00 %	\$	0.00 %

**Display method**

Print global figure only ☐  
List each client separately ☒

**Range Selection**

Scan all your employees ☐  
Scan only the employees from the selected State ☒

**Save** **Cancel**

Potential liability calculations are based on the Binaural percentage loss.

They can be run for a specific population of employees and return the total amount that would have to be paid in compensation if all the employees made a successful claim at the same time.

This figure is, of course, unrealistic but it is a good approximation of maximum risk.

If you have not entered each employee's state of origin, make sure you select "Scan all your employees"


For New Zealand, since there is no individual state, always have "Scan all your employees" selected.

Compensation figures are often indexed and change every year. Please make sure you obtain updated figure before running the report.

## Potential Liability

Make sure to enter the current figures for your state.

If your company operates in several states, you will have to run the report multiple times, or use the same figures for all..



**Australian Standard Setup**

**1- Not enabled**  
You can enter the audio Test data, and all reports use the thresholds values as entered.

**2- Enable AUSTRALIAN STANDARD - AS/NZS 1269.4:2014**  
You can enter the audio Test (first Test) plus a re-Test. Both can be done manually or automatically or mixed.

**3- Enable Modified Standard - (ESSO STANDARDS)**  
This is a variance on Australian Standard 1269. It shows the BEST thresholds instead of the MEAN at 3, 4 & 6 kHz.

**Australian Standard 1269.4:2014**

Enabled ☒  
Disabled ☐  
Enable Esso Standard ☐

**Display**

Show Retest figures on Audiogram report ☒  
Show Retest figures on Report card ☒  
Show 2000 Hz threshold as well ☐

**Effective Results**

Use Retest values at 3, 4 & 6 kHz for:


- % loss, Yes ☒ No ☐
- NALSO calculations, Yes ☒ No ☐
- Graphs on reports, Yes ☒ No ☐

**Page Setup** **Save** **Cancel**

## Australian Standard

If you select the Australian standard you will be able to store a second set of thresholds per test.

If you do, you may want to get the results printed on the Audiogram report.



**Browse Comments**

**Stored Comments**

Title	Comment
Abnormal test	Analysis of the above test show some abnormalities and we recommend that further tests be conducted to assess the matter further.
High frequency loss	High frequency sensorineural hearing loss both ears consistent with noise exposure. Binaural percentage loss of %. Percentage loss is calculated according to latest NAL formula. Low frequency loss is not considered work related.
Normal test	Analysis of the above test results show no abnormalities. A new test should be conducted in 2 year at the same time to conform to regulations.

**Add** **Edit** **Del / Undel** **Exit**

## Recording Stored Comments

You can record an unlimited amount of prepared comments ready to use, saving time and effort.

This will also have to effect of standardising your comments.

You can also store part comments and use cut and paste to build them up from pieces.

The software will allow you to compare the Class of noise Employees are exposed to with the class of the hearing protection device they are wearing.

The Class of noise the employee is exposed to is stored with the Workgroup. The Hearing protection device comparison mechanism in this software is based on The Australian Standard 1269.

**Warning:** This system is covers only Non Tonal noise levels up to 110 dB. For higher noise levels or tonal noises, please refer to a noise professional.  
The Protector Class will be automatically calculated based on the SLC80 value.

## 11.8 Setting up the Employee

---

Many different setting can be altered so that the software operates your way.



### Screen Layout

You can change some of the History questions. These questions exist both in **Employee** and **Test** in order to capture Medical history changes.



### Status, Tag & Re-train

Tags are useful to mark an Employee for some particular reason of your own. All the tags can be named for your own need. You can then do lists, reports etc. by tag. Tags can be removed globally as well.



### Special Labels

Click on the **Special Labels** tab and make sure the label names are suitable for your use. Each label gives meaning to a field which contains your employee's data. So if the label says **Employee Number**, that is what is meant to appear in the field, you may wish to change Employee Number to Bundy number for example. Click on **Save**.

## 11.9 Setting up the Testers

---

Goto **Setup >testers**. Add a tester and enter the name and if required the Approval number. You can add as many testers as required. Highlight the one most often used and click on **Set Default**. This name will be automatically loaded when adding a test.

## 11.10 Setting up the Test

---



### Test Setup

1- Select the **Percentage loss calculation formula** for your state. Normally you would calculate the percentage loss and correct for age (Presbycusis).

The program can use 2 different sets of tables for the calculation of the percentage loss. The Old table (1973) or the New tables (1988). The method of calculation will be indicated on the report.

In NSW compensation is currently calculated using the OLD tables. It is done using the NEW tables everywhere else.

If you do not want the Percentage Loss to print on the reports, this can be turned off in the **Audiogram** tab. See below.

2- Set the re-test gap for each type of Employee.

3- Alter the **Print out** fields to represent your company. Try a couple of printouts and alter the titles until you are happy with the result.



### **Test Questions**

1- In the **Test history**, a number of questions are identical to the one asked in the **Employee's history**, these are altered in the **Employee setup**. Some questions are specific to the test, they can be altered here.

2- You can customise the types of tests you will do. Type 4,5 and 6 can be altered.



### **Audiogram**

#### **1- Graph Colours**

The graph normal area and lines can be presented in any colour. If you have a black and white printer, select colours that will show properly.

#### **2- Audiogram Report**

You can turn On or Off any of the options shown in that box. This will affect the Audiogram report only.

#### **3- Threshold Normals**

The program can display a shaded area on the Audiogram marking the area of normality. The figures used to draw that line are also compared against the test results to produce an assessment of an Employee's hearing, and decide if the test is normal or abnormal. The Normal figures supplied with the program have been suggested by NAL and can be altered to suit your needs. If no normal value is required, enter 95 dB. The values can be altered in **Setup > Hearing Tests > Audiogram**

#### **4- Allow use of Mouse to draw audiogram**

To avoid errors by untrained staff, this feature can be turned off, specially if the data is transferred from an audiometer.



### **Spiro / Tymp / Bone**

In the **Spiro / Tymp / Bone** tab, un-check the option you do not need. If you select the Bone conduction, you will be able to import Bone conduction data directly from an Oscilla 960 or AD27 audiometer. Click on **Save**.

---

## **11.11 Setting up the Audiometer**

Go to **Setup > Audiometer**. Select the audiometer type connected to your computer. If you requested it, the software would have been supplied with a cable for your audiometer providing it is supported by the software.

Select the communication port, normally **Com1**, if using a USB to Serial converter, it will most likely by **COM 5**.

Add an Audiometer and enter Name, serial number and date calibrated. Add more audiometers if required. If this audiometer is the default, highlight it and click on **Set Default**. This name will be automatically loaded when adding a test. See also the Chapter on Audiometers.

Some audiometer come with their own software, or Noah software module. (Affinity, AS608, AD629e, MedRx, Astera, Aurical) In this case the program read files exported by the audiometer supporting software.

### 11.12 Australian Standard

---

Australian Standard 1269 requires that a test be conducted, followed immediately by a re-test at 3000, 4000 and 6000Hz . The real result would then be the mean of the two at 3000, 4000 and 6000Hz and the result of the test at the other frequencies. (see Australian Standard 1269 for exact details).

The software can store both these tests and then calculates the mean.

- 1- Australian Standard: the result is the MEAN of both tests
- 2- Esso Standard: the result is the BEST of both tests

There is no mention in the Australian Standard of how to calculate the percentage loss, so it is still calculated using the first test, but reports show the mean or best of both tests.

### 11.13 Potential Liability

---

The Potential Liability report is designed to give an estimate of the amount of compensation that would have to be paid in compensation, if all your employees made a successful claim at the same time. This figure is, of course, unrealistic but it is a good approximation of maximum risk. It could be used for financial planning, especially in case your company is self insured.

The report can be run for a specific population of employees. It can be run in two different ways.

- 1- using the relevant figures for your state, for all your employees, regardless of which state they are from.
- 2- using the relevant figures for one state, and scanning only the Employees from that state.

There are two formats:

- 1- report a global figure.
- 2- list each Employee and the compensation figure for each.

Because of the confidentiality of this report it is protected by a Master password, if you have activated the password function.

**Note:** the figures supplied are only for reference, you are responsible for entering the correct ones as they change normally every year.

### 11.14 Comments

---

The software can store any number of pre-defined comments. These comments can be loaded in the Employee's test at the time of the test with the click of a button. Once loaded, the comment can be further modified to suit your need.

The title of the comment is there just for your reference. When the **Comment** is selected only the first line of the comment shows, so the title is useful, to make a decision as to which comment to use. The title is not transferred to the Employee's test.

## 11.15 Letters

---

The software has a system of letters. You can create as many letters as you wish. They are saved in RTF format (Rich Text) and can also be edited by MS Word.

The software supports a simple but powerful form of Mail Merge. Fields can be selected when the letter is created. They are replaced with the relevant employee data at the time of printing the letter.

See the Chapter on Setting up Letters for more information.

Do not name your letters LAST\_LETTER.RTF, as it would be over-written.

## 11.16 Page Setup

---



### Report

You can set the top margins, the font, the number printed each time and the printer bin.

**Beware** that if you set the margins above 20, the report may not print when any component of the report is outside of the printable area. You would get a blank report. The font is common to letters and statistics



### Lists

You can select margins and font. Make sure you select a font that will produce an acceptable result with your printer, not all fonts are suitable.



### Letters

Margins are fixed here. To change the margins, use the **Margin Markers** in the letter.



### Statistics

You can change the margins and the colour of the Normal/Yes and Abnormal/No lines.



### Labels

This software is designed to use A4 sheets of sticky labels as available from Unistat, Avery and some other manufacturers.

Four pre-defined formats are available. You can also design your own by selecting the Custom label and entering the relevant figures in the label data window. The figures needed are :

- 1- the number of columns (1 to 3)
- 2- the width of one label
- 3- the height of one label

Labels print from top to bottom and from left to right.

**Note:** Continuous stationary labels are not supported as such, you will have to



experiment to see if you can design a custom label to match them.



### **MailMerge**

The mailmerge setup allows fields from the Client data to be selected for export to the **Dataset**. Click on the field to highlight it, then click on > to transfer the field to the Selected window on the right. In the same way fields that are not needed anymore can be removed by transferring them back to the left window.

If using MS Word, make sure that the **Field names needed** box is checked.

## **11.17 Printer**

---

This function allows the selection of any one of your printers. But in fact, only the default printer will be printed to. So to change the printer you will need to go to the **Control Panel** and make the new printer the default.

## **11.18 Preferences**

---



### **Backup drive**

You can select the drive where your backup will be saved to. This does not allow the selection of a specific directory. Data will be backed up to the root of that drive.



### **Advanced**

In some cases, for example if you have a network or a zip drive, you may want the backup files to end up in the sub-directory of a drive, not in the root directory.

This can be particularly convenient on a network when you want to do a backup every day in a different sub-folders called for example : F:\MONDAY, F:\TUESDAY etc ....



### **Running Another Program**

This feature is designed to allow running of another software or software module from within this software. To activate this function, go to **Setup, Preference,**

**Under Run another program 1 (or 2) ,** enter the name of the software and the title to appear in the **Files** menu.

for example:

Command line: C:\DUMMY\TEST.EXE

Title

Testing Software

Of course this software must be present in the correct directory before this function can operate properly.



### **Report**

You can define how many copies of the Audiogram will print each time and request that the all reports print to Screen first. You can then click on the **Print**, or **Zoom** button.



## CHAPTER 12 Setting Up the Audiometer

### 12.1 General Procedure

---

- Start the software and go to **Setup > Audiometer**
- Select the Audiometer type
- If the audiometer uses the Serial port, select serial port COM 1
- Add an Audiometer
- Enter the audiometer serial number and date calibrated, click on **Set Default** so that you do not have to enter the audiometer details for every test.
- Go back to the Main menu, Press the **Data entry** button
- Click on the **Employee** tab, Click **Add** and enter an employee, (fields with a red star are compulsory). Click on **Save**
- Click on the **Test** tab, click **Add**

#### **Serial Port Connected Audiometers**

- Click on the **Capture Test** button
- Initiate transfer (see each audiometer for particular)
- If transfer does not work, go back to **Setup > Audiometer** in the Software and change the serial COM port to another number.

#### **USB Port Connected Audiometers**

These audiometers are controlled by their own software, for example Calisto, Affinity, Medrex and others. Once the test has been completed by that software. The data can be exported and then captured by the Workplace Audiometrics software.

- Start the original Audiometer software. Enter a test and export the test right away (before saving it) .
- Click on the **Capture Test** button
- Select the file you just exported

### 12.2 My Laptop does not have a Serial Port

---

If your audiometer is one of those above who needs a serial port, but your computer does not have one, please follow these additional steps:

- Purchase a USB to Serial converter (we do carry these)
- Install the USB to Serial driver from the CD that came with it. When requested to do so...
- Connect it to a USB port on your Laptop, connect the output of the converter to the cable provided with the software and connect the other end to the audiometer.
- Go to **Setup > Audiometer** and select the serial port COM 5 instead of COM 1 and try again as above. Or go to the **Control Panel > System > Device Manager** and find out what Serial port your converter is set to. Then select that number in **Setup > Audiometer**.

### 12.3 Installing a USB to Serial adapter

---

If your laptop does not have a serial port (most recent ones do not) you will have to purchase a USB to serial adaptor. They are inexpensive and connect into a USB port on your laptop.

This adapter will only work after you install the driver software. Beware that the driver software comes on a small mini CD **hidden inside the packaging**. Make sure to remember the location of the USB port you use as the driver will only work with that port.

Once the driver is installed on to the laptop go to the Control panel, select **System > Device Manager > Universal Serial Adaptor**. Find your adaptor and check which serial port it is connector to.



If needed you can change the port number. This port number will be required to be entered in the Workplace Audiometrics program under **Setup > Audiometer > Com Port**.

Remember that you still need to use the gray interface cable supplied with the software from the adaptor, to the audiometer.

Each audiometer uses a differently configured cable, so it is important to use the cable designed for the audiometer you are using.

## 12.4 The OSCILLA SM910 / 930 / 950 and 960

These audiometers are very simple to use, after entering the Employee information, add a new test into the computer, click on the **Capture Test** button (or **Capture Retest**) in the software. This will put the computer in waiting mode.

- Press the **AUTO** button on the audiometer or perform the test manually.
- When ready, press **DATA** on the audiometer to send the data to the waiting software on the computer.

OSCILLA			
COMPUTER			AUDIOMETER
female	DB9		DB9 female
TXD	3	-----	2
RXD	2	-----	3
GND	5	-----	5
DCD	1	-----	1

## 12.5 The Maico MA30

**Setting up:** Connect the male end of the cable provided to the serial port at the back of the MA30. Connect the other end of the cable to the serial port at the back of your computer.

- On the MA30, Press the **Menu** button
  - Move the cursor down to "**PC or Printer**", press **Enter**
  - Move up to "**PC Data Transfer**", press enter
- Press **Menu** to exit the menu system.

This audiometer is very simple to use, after entering the Employee information, add a new test into the computer, press the **Capture Test** button.

- On the audiometer perform the test either manually or automatically. The results will be

displayed in the software window. When the test is complete, press on the **Save** button.

MA30			
COMPUTER		AUDIOMETER	
female	DB9	DB9	male
RXD	2	-----	2
RTX	3	-----	3
GND	5	-----	5
DSR	6	-----	6
CTS	8	-----	8

## 12.6 The Maico MA15

---

- The MA15 must be set to **print to the Computer**, not to Printer. Setup is entered by pressing Help while turning power on. Use patient button to scroll down.
- You must have the original printer cable from Maico
- Connect the other end of the MA15 printer cable coming from the audiometer to the cable supplied with the software. Connect the other end of the cable supplier with the software to your computer's serial port.

To transfer the data, wait for the end of the test, then press **Print** twice.

MA15			
COMPUTER		AUDIOMETER	
female	DB9 DB25	DB9	female
RXD	2 -- 3	-----	1
RTS	7 -- 4	-----	2
GND	5 -- 7	-----	3
[	4 20		
[	6 6		

## 12.7 The InterAcoustics AD27

---

The AD27 can be used in manual transfer mode (one threshold at a time or in automatic screening mode. The option is selected by selecting the correct audiometer from the list ie AD27 manual

### Manual Transfer:

- Press Store on the audiometer to transfer the current frequency and threshold.

### Automatic test transfer:

- Start the Audio test on the audiometer by pressing "Audio Threshold"
- When the test is completed, come back to the computer and click on the **“Capture Test”** or **“Capture Re-test”** button.

AD 27			
COMPUTER		AUDIOMETER	
female	DB9 DB25	DB25	male
TXD	3 -- 2	-----	2
RXD	2 -- 3	-----	3
CTS	8 -- 5	-----	8
GND	5 -- 7	-----	5
DTR	4 -- 20	-----	20
		-----	6

## 12.8 The InterAcoustics AS 216 and AD 226

---

To setup the audiometer press **Shift** while turning the **Frequency** knob.

**Baud Rate:** Turn the knob until # 1. Turn the “**HL dB**” knob to read 38.4 (38400 Band) lower Baud rate are ok too but must be reflected in the software.

**Hardware handshake:** Turn the Frequency knob until # 12. Turn the “**HL dB**” knob to read “**HAOF**”

**Transmit on Store:** Turn the Frequency knob until # 13. Turn the “**HL dB**” knob to read “**trOF**”.

**Close Setup:** Press Shift while turning the Frequency knob. Setup is complete.

To capture test data:

- Perform the test in Manual or automatic mode. When the test is completed
- In Data Entry, Add (or Edit) Test, click on the “**Capture Test**” or “**Capture Re-test**” button, then click on the “**Start**” button. This will capture the data including Bone if Bone has been tested.

AS 216/ AD226				
COMPUTER			AUDIOMETER	
female	DB9	DB25		DB9 male
TXD	3	2	-----	3
RXD	2	3	-----	2
CTS	8	5	-----	8
GND	5	7	-----	5
DTR	4	20	-----	4

## 12.9 The GSI 66 / 67 / 68 Audiometers

---

These audiometers must be setup properly in a specific manner. Further details can be found in the audiometer manual under Setup. In brief to access the Setup mode: turn the power on while pressing the X button.

- Press Enter 5 times
- Memory: Press X to toggle between ON and OFF, press Enter
- Setup, press X until you see **Printer Computer**, press Enter
- Press Enter
- Level, turn Tone knob until you get -5 dB, press Enter
- Frequency      125 , press X until it is OFF  
                          250, OFF  
                          500, ON  
                          1000, ON  
                          1500, ON  
                          All others ON, press Enter
- Tone , Press Enter
- H factor, press Enter, press Enter to exit Setup.

The data is transferred to the computer at the end of the test by pressing Print then Enter

GSI 66/67/68				
COMPUTER			AUDIOMETER	
female	DB9	DB25		DB25 male
TXD	3	2	-----	2
RXD	2	3	-----	3
RTS	7	4	-----	4
CTS	8	5	-----	5
GND	5	7	-----	7

## 12.10 Audiometer Data Transfer

---

The software presently supports quite a few audiometer types. The Microlab, MA800 and MA728m have a special feature. They can store tests internally and dump the entire storage into

the software.

All the above Audiometers can be used in direct mode:

- The Employee is present at the test site
- The Employee data is entered into the program while the test is being done
- The data is transferred from the Audiometer to the Computer at the end of the test.

The Microlab, MA800 and MA728 can be used in memory mode:

Part 1

- You are in the field without a computer
- The tests are done and recorded in the audiometer's memory
- You record on paper the Audiometer sequence number and the Employee details

Part 2

- Upon returning to base, you connect the audiometer to the computer
- You tell the computer which Audiometer you are using  
Go to **Setup, Audiometer** select the correct type
- Transfer all the tests in one go to the computer

Go to **Files, Download from MICROLAB** or **Download from MA800/728**

Part 3

- One by one you identify and find (or enter) the Employee's details using the **Capture from File** button. Use the sequence number as a match, each time, reading the Employee test from the audiometer file.

---

## 12.11 Reading from the Earscan or Microlab

---


The transfer can only work with the audiometer command to send data to the printer like GRAPH, DATA on the EARSCAN and 7, 0 7 Y on the MICROLAB. The commands to send data to a computer like S (MICROLAB), will not work and only produce garbage.

---

## 12.12 Earscan / Microlab Cable Diagram

---

A D25 or D9 socket is used depending on the type of connector on the computer. The cable can be maximum 10 metres long.



---

## 12.13 Using an Earscan with Air and Tymp

---

Press IMP to do the Impedance test, then press AUD twice to do the air test.

When the test is completed:

- Select the Employee
- Add a test, enter the Test details in the computer.

When the program requests the method of acquisition of Audio Data:

Press DATA: Only the figures will be sent (no curves).

The output of the I=TYMP command is exactly what would come out of a serial printer connected directly to the Audiometer. Each time a new graph is transferred the old graph is over-written. Only one graph is ever present in the computer.

## 12.14 Using an Earscan or Microlab with Air only

---

Press AUD twice to start the automatic test. When the test is completed, select the Employee, add a test, enter the Test details in the computer. When the program requests the Method of Acquisition of Audio Data, select EARSCAN then press DATA on the audiometer.

## 12.15 Single test transfer with the Microlab

---

Select the Employee, add a test, enter the Test details in the computer.

When the program requests the method of acquisition of audio Data, select **Microlab**.

On the **Microlab**:

- Press **5- Recall Test**
- Press **3- Request ID number**, enter the ID number
- Press **7- Transfer**

...and so on for each test. In this way you have control of the order of the transfer.

## 12.16 Multiple test transfer with the Microlab

---

The Microlab has a memory buffer to store up to 150 tests. When you want to transfer the stored tests out of the audiometer, two methods are available:

The entire data stored in the **Microlab** can be dumped into a file, and the file used to supply the results to the Employees's test one at a time in sequence.

- Go to the **FILE** menu
- Select **Download from Microlab**

On the **Microlab**

- Press **0**
- Press **7**
- Press **Y**

Now when you want to enter tests, the "Capture Test" and "Capture Re-Test" buttons have changes to "Test from File" and "ReTest from File".

The idea is to read the results from the file downloaded from the Microlab, one at a time. A label will display the last test number read. Hopefully you would have recorded the name and test number when doing the tests, so now, you enter or go to the existing Employee's record, add a test and read the thresholds from the file in the corresponding order.

- Go to the data entry screen
- Find (or enter) the Employee corresponding to the test
- Add a test, press the "Test from File" button to extract the data from the file into the test.

The next ID number will now show. When the file is empty of tests, the ID number will disappear and the buttons will revert to "Capture Test" & "Capture Re-test".

To cancel the “Test from File” mode, go to the File menu, Select Download from Microlab and Exit without sending any data.

---

### 12.17 Multiple Test Transfer with the MA800 / MA728

---

The entire data stored in the MA800 OR MA728 can be dumped into a file, and the file used to supply the results to the Employees's test one at a time in sequence.

- Go to the FILE menu
- Select Download from MA728

On the Audiometer

- | MA800                 | MA728                       |
|-----------------------|-----------------------------|
| - Turn Audiometer ON  | - Turn Audiometer ON        |
| - select TESTS & BL'S | - Press 8 (Data transfer)   |
| - select TESTS        | - Press 61 (Transmit Tests) |
| - select XMIT         |                             |

The data will be dumped in a file called AS\_XMITM.DAT

Now when you want to enter tests, the “Capture Test” and “Capture Re-Test” buttons have changes to “Test from File” and “ReTest from File”.

The idea is to read the results from the file downloaded from the Microlab, one at a time. A label will display the last test number read. Hopefully you would have recorded the name and test number when doing the tests, so now, you enter or go to the existing Employee's record, add a test and read the thresholds from the file in the corresponding order.

- Go to the data entry screen
- Find (or enter) the Employee corresponding to the test
- Add a test, press the “Test from File” button to extract the data from the file into the test.

The next ID number will now show. When the file is empty of tests, the ID number will disappear and the buttons will revert to “Capture Test” & “Capture Re-test”.

To cancel the “Test from File” mode, go to the File menu, Select Download from Microlab and Exit without sending any data.

---

### 12.18 Setting up the Microlab / MA728 / MA800 to store data automatically

---

These audiometers can store multiple test for later recall.



#### **Microlab**

The storage must be requested for each test by pressing the STORE TEST key.





## MA800 / MA728

### MA800

- Turn the audiometer OFF then ON
- Select SETUP
- Select STORAGE
- Select AUTOMATIC TEST STORE ON
- Exit

### MA728

- Turn the audiometer OFF then ON
- Press 6 (Storage Options)
- Press 21 (Auto test storage)

## 12.19 Setting up the Maico MA800

- Turn the audiometer ON
- Select SETUP
- Select DATA TRANSFER and set the following
 

BAUD RATE	9600	HARDWARE HANDSHAKE	ON
PARITY	NONE	DELIMITED DATA	ON
STOP BITS	1	EXTENDED DATA	OFF
XON/XOFF	OFF		
- Press EXIT and turn audiometer OFF
- Turn audiometer ON, it is now ready to be controlled by the computer. During the test the entire operation is controlled by the computer.

### MA800

	COMPUTER		AUDIOMETER
	female DB9	DB25	DB25 male
TXD	3	2	2
RXD	2	3	3
RTS	7	4	4
CTS	8	5	5
GND	5	7	7


In some cases to get the new MA 800 to work we had to select Even parity on the audiometer. Bear that in mind if you experience communication problems like garbage coming out.

## 12.20 Setting up the Maico MA728m

This audiometer operates in the same way as the MA800 but the cable and the communications parameters are different.

The MA728 MUST be setup internally to operate at 9600 Baud. It normally comes set at the factory to 1200 Baud. To alter the settings, the Audiometer cover must be removed, the pc boards cover plate removed by pinching open the 5 black clips and lifting them.

The communication board is the third from the left (see diagram below)



Do not remove the communication board, look at it from the left of the audiometer. It has 2 set of vertical switch rows. Do not alter the setting of S2. The top switch of S1 must be moved to the right and all others on that vertical row must be to the left. When this is completed, reassemble the audiometer.

MA728				
	COMPUTER			AUDIOMETER
female	DB9	DB25		DB25 male
	3	-- 2	-----	2
	2	-- 3	-----	3
	7	-- 4	-----	4
	8	-- 5	-----	5
	5	-- 7	-----	7
	4	-- 20	-----	6
	1	-- 8	-----	8
	6	-- 6	-----	20

## 12.21 The Rexton AS53, AS54 & DA63

---

This audiometer operates in a similar fashion as the Earscan audiometer. A similar message appears on the screen when transfer is required. All is needed is to select the correct audiometer from the list in the SETUP menu, AUDIOMETER option.

- Enter the Employee details, add a test and enter the relevant details, when asked the mode of entry of the Audio details, select **AS53/54**. The computer will then be waiting for the audiometer to send the data across.

- Perform the test either automatically or manually. (by pressing the **MAN/AUT** button.)  
When the test is completed, press the **X-mit** button to transfer the data to the waiting computer.

AS 53 / 54, DA 63				
	COMPUTER			AUDIOMETER
female	DB9	DB25		DB9 male
TXD	3	-- 2	-----	3
RXD	2	-- 3	-----	2
RTS	7	-- 4	-----	7
CTS	8	-- 5	-----	8
DSR	6	-- 6	-----	6
GND	5	-- 7	-----	5
DTR	4	--20	-----	4

## 12.22 The Tracor RA 300 / 500

---

This audiometer is very simple to use. After entering the Employee information, add a new test into the computer, select from RA500, perform the test on the audiometer either in automatic or manual. When ready, press the “Capture Test” or “Capture Re-test” button.

RA 300 / 500				
	COMPUTER			AUDIOMETER
female	DB9	DB25		DB25 male
TXD	3	-- 2	-----	2
RXD	2	-- 3	-----	3
RTS	7	-- 4	-----	4
CTS	8	-- 5	-----	5
DSR	6	-- 6	-----	6
GND	5	-- 7	-----	7
DCD	1	-- 8	-----	8
DTR	4	-- 20	-----	2

## 12.23 The Tracor RA 400

As for RA 300

RA 400			
COMPUTER		AUDIOMETER	
female	DB9 DB25	DB25 male	
TXD	3 -- 2	2	
RXD	2 -- 3	3	
RTS	7 -- 4	4	
CTS	8 -- 5	5	
GND	5 -- 7	7	

## 12.24 The AudioScan 1000

The serial cable attaches to the Printer connector at the back of the audiometer.

To transfer the results press the **Capture Test** button in the software, then start the test either in manual or automatic mode.

After the test is completed, just press the **Print** button on the audiometer then press **Yes**.

If **Auto Print** is on, the transfer will happen automatically if the operator has pressed the **Capture Test** button prior to starting the test.

AS 1000			
COMPUTER		AUDIOMETER	
female	DB9 DB25	DB9 female	
RXD	2 -- 3	2	
TXD	3 -- 2	3	
DTR	4 -- 20	4	
GND	5 -- 7	5	

## 12.25 The GN Resound Xeta & Danplex 300a

**Setting up:** On the Audiometer press the **Setting 1** button so the light on the switch is on. Then press the **Setup** button.

- Rotate the **left dB knob** until the display indicates **Baudrate**.

- Rotate the **right dB knob** until the display indicates **57600**.

- Press the **Setup** button again to save the settings.


If you are also using **Setting 2**, repeat the operations with **Setting 2** button illuminated.

In the Workplace Audiometrics program, go to **Setup > Audiometer**, select the Xeta audiometer and tick the same Baudrate you selected earlier (57600) then select the com port (most likely Com 1). If you are using a USB to Serial converter, the Com is likely to be Com 4 or Com 5.

When the patient responds to a tone presentation, the corresponding threshold value is shown on screen in the correct Frequency column, and memorised. Air and Bone conduction thresholds are stored automatically, but masking information is not retained. The tests can be done manually or automatically.

## 12.26 The AS 608, AD629e, Affinity, Calisto & Astera Audiometers

These audiometers are controlled by their own software, sometimes even as a module under Noah. The test is done with that software, then imported into the Workplace Audiometrics software.



Procedure to Import a test:

Once the audiogram is complete, don't save, simply click **Menu > Edit > Export**


You will then be prompted to save the file.

- Locate a folder you wish to save the file to, name it and save with .xml as extension.
- The file is then ready to be captured by the Workplace Audiometrics software.
- Go to the Workplace Audiometrics software
- Add a test, click on the Capture button
- Select the file name created above

If you always capture the test right after a test, you can use the same filename every time.

Instructions for older saved audiograms:

The only problem we ever have with importing and exporting audiograms is that customers export previous sessions, in these cases the xml files are blank.



Open the audiogram you wish to export which has been saved.

Select **Menu > Edit > Transfer to Current Session**.

Then click **Menu > Edit > Export**

#### Setting up:

You can tell the Workplace Audiometrics software which audiometer you want to use and where to find the exported files.

Go to **Setup > Preference > Audiometers** and select the drive and folder there.

---

## 12.27 The MedRx Audiometers

The Medrx Avant A2D works like the audiometers just above

In the Audiogram screen of the MedRx program, select **File > Export**, then save the file with the extension .txt.

Make sure to use the **Export** option, not the **Save** option. The Workplace Audiometrics program is not able to read the files generated by the Save option, only the files generated by the Export option.

---

## 12.28 The Aurical Audiometers

- 1- Perform the test with the software
- 2- Click on Edit, then on Export
- 3- Type a file name, making sure it has the extension of .xml Example **Test.xml**
- 4- Go to the Workplace Audiometrics software
- 5- add a test, click on the Capture button
- 6- Select the file name created above

If you always capture the test right after a test, you can use the same filename every time.

You can tell the Workplace Audiometrics software where to find the exported files. Go to **Setup >Preference > Importing Data** and select the drive and folder there.

## 12.29 Trouble shooting Communication Problems

---

- ◆ The major cause of problems initially is the selection of the wrong serial port. This can be fixed by changing the port using the SETUP menu, SET AUDIOMETER option, or by connecting the audiometer to the other port.

You may have a mouse on COM1, try COM2. The serial port on your computer is a male plug, it may have 25 or 9 pins. Some computers have a 25 pins as COM1 and a 9 pins as COM2 or vice versa. COM3 and COM4 are a rare occurrence and they normally prevent the proper use of COM1 or COM2 because they share the same interrupt.

- ◆ The wrong parameters have been used either in the program or in the audiometer.
- ◆ With the EARSCAN one problem may be that the output is set for COMPUTER, when it should be set for PRINTER. When the Machine is first turned ON, a message appears:  
MAKE SELECTION OUTPUT TO PRINTER  
If the message is OUTPUT TO COMPUTER, then it is incorrect.
  - Press SPEC, then press DATAThe message should now be correct (OUTPUT TO PRINTER)
- ◆ With the MA800/MA728 a main power spike may have upset the synchronisation between the computer and the Audiometer. Exit from the data capture screen then turn the Audiometer OFF then ON again.



## CHAPTER 13 The letter System

### 13.1 Principal of Operation

---

The software includes a powerful, easy to use letter system. Letters are stored as **Form Letters** in RTF (Rich Text) format. The form letters are used to create the **Employee Letters**.

To create and print an Employee letter

- go to the Data entry screen and click on the **Print** button and select **Letters**
- select an appropriate letter from the list of RTF files

The **Fields** in the letter will be expanded with the content of the Employee's information and the letter printed. The letter will be saved as Last\_Letter.RTF, so that you may review it later. Only the last letter saved can be retrieved as the same name is used over and over again. The original Form Letter is not altered by this process.

### 13.2 Creating / Altering the Form Letters

---

Go to **Setup, Letters**, Type some text then click on the Letter menu and **Save As** or **Save**. You may also open an existing letter to use as the basis for a new letter, then **Save As** with a new name.

### 13.3 Inserting a Replaceable Field

---

Replaceable fields are inserted in the **Form Letter** when it is created or altered in the **Setup** menu.

- Move the cursor where you want to insert the field
- Click on the **Insert Field** button, a list of the fields available will appear
- Use the mouse to highlight the field you want then click on it to load it.

A name will be inserted as in **[Employee\_Name]** . This will be replaced by the real Employee's name at the time of printing. Example: Creating a full name:

- Select Name, Surname, from the field's list, one after the other.

example: **[Employee\_Name] [Employee\_Surname]** is inserted

result at print time: "Paul SMITH" (the Employee's name)

### 13.4 Creating Gender Sensitive Sentences

---

Use the replaceable fields

- Select His/Her

example: {[HIS/HER] problem is...}


result at print time "His problem is..." (if the Client is male.)

- Select Mr/Mrs, Name, Surname


example: **[MR/MRS] [Employee\_Name] [Employee\_Surname]**

result at print time: "Mr Paul SMITH"

The field's name is enclosed in special characters for recognition by the programme at the time of replacement. Should you delete one of those characters you can remove the entire field's name and re-select it from the list.




In the image above you can see the Merge Fields pull down menu revealing the list of fields available. More fields are listed which are not visible on this image.




You can view page breaks by pressing on the Thumbnail button.

There are no **Print Preview** options, because what you see on the screen is exactly what you will get on the printout.


The letter editor closely emulates the operations of MS Word.



The final letter above is based on the template from the previous page. It contains all the relevant Employee's data. The red squiggles are created by the spell checker to highlight possible spelling errors or data it does not understand.



To update old versions of the letter templates, if you have any, use the button pointed to above, then move the ruler markers as shown below, before saving the file.





## 13.5 Using Bold and Underline

---

You can use any of the font attributes available in you text, this include **Bold**, Underline, *Italic*, font size, font type and font colour. Font colour is selected from the **Edit** menu, **Font** option.

## 13.6 Changing the Letter Margins

---

The margins can be changed by moving the left and right margin arrows (Margin Markers) into the right position. If there is text present in the letter you should highlight the text first. The top left Margin marker is the left margin for the first line of a Paragraph.

**Standard window envelope format** : 230mmX100mm and the window is 58mm from the top and 32mm from the left side.

The default margins should be

Top	= 20
Bottom	= 20
Left	= 24
Right	= 24

This can be altered in **Setup > PageSetup > Letter**

The letter system allows the creation of template which are used to create Employee letters very quickly with the press of a button.

Creating a template: the operator, in the screen below has selected to insert a Merge field for the Employee's surname. This is done by selecting from a list of merge fields in the pull-down list.

You will notice other merge fields in the letter. They are enclosed in square brackets like **[Employee\_Name]**. In some cases there are both upper and lower case version of the merge field.

There are also pronouns to be used to make up correct sentences like "**His** hearing... or **Her** hearing.... As well as **his** and **her** for use in mid-sentence.

Once finished, the template is saved for later use. The template letters are used against an employee in the **Data Entry** screen, or against a group of Employees in the **List Employees and Tests** screen. All the merge fields are replaced with the correct employee's information.


In the **Data Entry** Screen, the text of the letter can be altered as required before printing or saving.

In the **List Employees and Test** screen, the letters will be printed as per the template.

## CHAPTER 14 The Spell Checker

The spell checker include two functions: Spell checking on demand (through the **Spell** button) and Live Spell checking, which is active by default.

### 14.1 Setting up the Spell Checker



Go to **Preference > Option** and click on the **Configure** button.

The default dictionary is **British.adm**, but more dictionaries can be added, there is one available for American English called **American.adm**. Other languages are available on request.

The Spelling checker offers a custom dictionary which is added to when you select **Add a word** from the Spelling checker dialogue.

The default custom dictionary is called **Spell.adu**, there is also a learning dictionary call **Learn\_sp.adl**


Live spelling can be de-activated by un-ticking the **Live Spell** checkbox.

### 14.2 Spell Checking a letter



S

In the letter, click on the **Spell** button. When a miss-spelt word is encountered, the pelling dialogue will appear automatically.



You can select a action from the buttons. To add the highlighted word to you own dictionary Click on the **Add** button. The dictionary can be changed using the **Options** button.

The spelling checker is available and active in all the memo fields and comments of the program as well as in the letters.

Live spelling is also active in these areas but can be de-activated by un-ticking the **Live Spell** checkbox.



## CHAPTER 15 Principles of Operation

### 15.1 Percentage Loss

---

The percentage loss is extracted from the values of the hearing levels of the subject at the audiometric frequencies 500, 1000, 1500, 2000, 3000, and 4000 Hz. The age is used to calculate the presbycusis correction (loss of hearing due to aging)

**1500Hz:**

If the audiometer is not fitted with 1500 Hz and BOTH the left and the right thresholds at 1500Hz are absent, a mean threshold is obtained from the value of 1000Hz and 2000Hz and used for the calculations.

If 1500Hz is only missing on one side the % loss for that side will not be calculated.

The program can use the 1973/1976 NAL table or the new 1988 formulas for the calculation of the % hearing loss. The binaural % loss with the old table is calculated using the following formula:  $(\text{best ear} \times 4/5) + (\text{worst ear} \times 1/5)$ .

**Missing frequencies:**

If one or more of the six frequencies used to calculate the percentage of hearing loss other than 1500Hz, has not been entered, the result will be replaced with '--.' as there is insufficient data to obtain a valid result.

**Threshold greater than maximum:**

If one or more of the hearing levels are greater than the maximum preset dB level (indicated by 'V') the corresponding percentage of hearing loss will be preceded by the sign 'greater than' EG: >50.0%

**Compensation in NSW:**

Normally the 1988 NAL tables are used. In NSW the law used to requires that compensation % loss be calculated using the 1973/1976 tables with the compensation method of calculating the Presbycusis (.5 dB per year above 50 years old).

**New ruling for NSW:** From January 2002 Workcover has issued a general recommendation to use the 1988 table for compensation in NSW.

The Option menu in the SETUP allows changing of the Tables and methods of calculating the Presbycusis. The % loss can also be withheld from the report if required.

### 15.2 Presbycusis Correction

---

The PRESBYACUSIS is the normal amount of hearing loss resulting from the aging process. After the total % loss is calculated, the presbycusis component is subtracted from the total % loss to obtain the corrected value. This process is applied to the left, right and binaural values.

The Presbycusis can be calculated using the corresponding tables/formulas (NORMAL) or using the NSW Compensation method (COMPENSATION). This method is used in conjunction with the 1973 tables for compensation cases in NSW.

In the Compensation method a correction of .5dB per year above 50 years of age is applied to each threshold before calculation.

Example:  
Frequency 500 1000 1500Hz Age=70  
Threshold 35 45 55dB | correction is  
corrected 25 35 45 ← (70-50) x .5 = 10  
% loss 1976 1.5% +5.2% +7.4%

In the Compensation calculation the binaural loss is calculated using:  
(best ear x 4/5) + (worst ear x 1/5).

Example:      LEFT   RIGHT      BINAURAL  
                 20%   15%      16%

Best ear is Right ear →  $(15 \times 4 / 5) + (20 / 5) = 16\%$

### 15.3 The % Loss Report

---

A special report allows for presentation of a Percentage Loss chart for compensation purpose in Australia.

### 15.4 Normals Assessment

---

To establish if a subject's hearing is normal, the results are compared with a set of normal values. These values can be altered in the **Setup > Hearing Test > Audiogram** menu. The default values have been suggested by NAL.

If the limit is greater than the maximum available on your audiometer (95dB) then the limit is cancelled.

The normal values are not corrected for age as no standard exists for this.

**Note:** Because the threshold data can only be entered in 5dB steps, a borderline normal (but still in the normal value range) might get up to 3% percentage loss. The normal limit should really be set to 18dB. So do not be surprised if a "normal" person still has a small loss.




## CHAPTER 16 Interpreting the Audio Test Results

### 16.1 Assessment of the Audiogram

Audiometric results can be a little confusing to interpret. Below are some examples of confusing reports and how to understand them. These reports by way of example, may shed light on the subject.

#### The Audiogram Report

Consider the following Audiogram Report



The questions regarding the above report were:

**Q-** Why is the percentage loss on the left = 0, when the thresholds at 6kHz and 8kHz are so elevated?

**A-** The percentage loss (except in Western Australia) does not take into account the thresholds at 6 and 8kHz. If you are testing in Western Australia, you need to change the formula used for the calculation. You will find this in Setup > Hearing Test > Test Setup.

For Western Australian select Formula 3, for other Australian states set it to Formula 2.  
 Formula 1 is no longer used in Australia since 2010.

**Q-** Why is there a percentage loss on the right, when the thresholds are normal?

**A-** The normal values supplied by NAL and stored into the program are used to assess if a test is normal. You can change these values in Setup > Hearing Test > Audiogram.

At 1kHz, the normal is 20dB, however the percentage loss at 18dB at the same frequency indicates a loss of 0.2%.

The problem is that audiometers test thresholds in 5dB steps. So you can test 15dB or 20dB, but not 18dB. So 20dB which is very close to 18dB indicates a small loss.

**Q-** Why is the normal level for 8kHz set to 95 db?

**A-** Hearing sounds at 8Khz have little benefits in everyday work life, unless you are a musician. So the level is set to 95dB so as not make the test abnormal and to raise the alarm.

**Q-** Why does the Australian Standard bother to test 8kHz then?

**A-** A loss at 6 or 8kHz may well indicate that the individual has been exposed to loud noises in the past, and as a result has suffered a permanent loss at the high frequency end of the spectrum.

**Q-** Where does the normal values come from ? (Shaded area on the audiogram)

**A-** The assessment of normality of the test is made based on normal values provided by National Acoustics Laboratory. These defaults can be altered in the software in **Setup > Hearing Test > Audiogram** if they do not suit the user.

Date	Left ear								Right ear									% Loss		
	0.5k	1k	1.5k	2k	3k	4k	6k	8k	0.5k	1k	1.5k	2k	3k	4k	6k	8k		Left	Right	Binaural
07/09/2000	10	0	--	10	25	35	35	25	5	0	--	5	15	30	20	40	B	2.2	0.8	1.3
	Re-test -----> -- -- --								-- -- --											
11/07/2001	15	+5	--	10	20	30	30	35	10	+5	--	+5	10	25	25	45	Y	1.1	0.3	0.7
	Re-test -----> 25 35 25								15 35 30											
16/08/2002	15	0	--	+5	25	35	25	25	10	+5	--	+5	10	25	25	45	Y	2.2	0.3	1.0
15/11/2010	15	10	10	10	35	45	40	40	20	10	15	20	35	50	45	40	B	5.2	7.2	6.0
	Re-test -----> -- -- --								-- -- --											
10/09/2013	15	5	10	10	30	40	60	40	15	5	5	10	25	40	65	55	Y	2.6	1.9	2.2
	Re-test -----> 30 40 60								25 50 65											
12/09/2013	10	5	5	5	25	45	40	45	10	0	5	10	25	40	50	50	R	2.9	1.9	2.1
	Re-test -----> 25 45 40								25 40 50											
	0	0	0	0	0	0	0	5	0	0	0	0	0	0	5	10	Loss (last test - baseline)			

Number of tests: 14

**Australian Standards**

(a) a shift in mean threshold at 3, 4 & 6 KHz >= 5 dB  
 (b) a shift in mean threshold at 3 & 4 KHz >= 10 dB  
 (c) a shift in thresholds at 6 KHz >= 15 dB  
 (d) a shift in mean threshold at 0.5, 1, 1.5 & 2 KHz >= 15dB  
 (e) a shift in thresholds at 8 KHz >= 20 dB

(16,1) Victorian Regulation 2007 3.2.12 & NZ Std  
 (f) current mean thresholds at 1, 2 and 3kHz => 25 dB (Rio)  
 (h) a shift from baseline av. At 2, 3 & 4kHz =>10dB (Rio)  
 (1) a shift from baseline av. at 2k, 3k & 4kHz => 10 dB, age adj.  
 (2) Average abs. at 2k, 3k and 4k => 25 dB same ear as in (1)

Call center operator: Fail

Train driver: Pass

(a) (b) (c) (d) (e)

☒ ☒ ☒ ☐ ☐

(16,1) (f) (h) (1) (2)

☐ ☐ ☐ ☐ ☐

NIHL last test

120

125

## 16.2 Assessment of the Report Card

Consider the following Report card

**The questions regarding the above Report card usually are:**

**Q-** Why does the Australian standards indicate a shift for (a), (b) and (c) when the loss on the last line only indicates 5 db at 6kHz and 10 at 8kHz?

**A-** The last line only indicates a loss of hearing (Hearing has worsened), the Australian Standard only mentions a shift. In this case it is a positive shift (hearing improved)

You will notice that there are 2 baseline tests, one in 2000 and one in 2010, only the last one

(2010) is taken into consideration.

**What is the best report for detecting the onset of hearing loss?**

**Q-** What report is the most helpful in detecting the onset of hearing loss?

**A-** There are 2 report which indicate a change in threshold in groups of people

**1- The baseline report**

This report will scan group of people and flag those whose have a loss  $\geq 15$ dB at any frequency between 2 and 6kHz

This is a good report, however the time span between the baseline and the latest test may be quite large. Therefore this might trigger a false alert, for example of the baseline test was tested in 1990 and the last test was tested in 2013, in 23 years 15dB hearing loss may be acceptable depending on the age of the person.

**2- The Reg.2007.3.2.11 report**

This report is compulsory in Victoria, however it is a great report as it will flag a loss of 15dB since the last test. There may be one or two years between these and a 15dB loss in 2 years should act as a warning and trigger some vigilance or retest to assess if there is a problem.





## CHAPTER 17 - Network Use

To use the program on a network you must have the Network version of the program and a site licence to cover a number of additional Workstations. A network is composed of a Server computer connected to one or more Workstation computers. If you only have the single user version, you can also upgrade to the network version. You will notice that the main screen now shows : **Network node = xxxxxxxx** at the bottom right of the main screen.

### 17.1 Installation

---

The software is installed on one host computer only, on the network. This computer must be visible to all the other computers on the network, and the folder where the program is installed must be shared on the network. In addition all users must have read/write access.

### 17.2 Host Setup

---

- Once the software is installed on the host,
- using Windows explorer right click on the program folder (C:\Audio)
  - Click on Properties > Sharing > Advanced Sharing
  - Tick Share this folder
  - Click on Permissions and give all Full control

### 17.3 Other computers Setup

---

There is no need to install the software on the other computers. It is only installed once on the host computer.

- using Windows select **Tools > Map Network Drive** or **Computer > Map Network drive** (Windows 8)
- For Drive, select the K: drive (for example)
- Click on the Browse button and select the folder you shared earlier on the host
- Tick Reconnect at logon
- Finish

Adding a shortcut on the other computer

- Click on the desktop, select - New > Shortcut

Command line : Enter the network drive followed by the directory then the program name

Example K:\SCREENW.EXE

Next select a name for the short cut?

- SCREEN
- Finish OK

When you click on the icon the program will be loaded from your host and run on the local computer.

### 17.4 Parameter Change

---

The parameters you setup are common to all the computers on the network.

When making changes to a record or adding a record, the program locks the record or file temporarily preventing other users from accessing the data. When this happens, the program

retries for 5 seconds before asking if you want to try again as the other user may have finished with the record by then.

The parameters are loaded by each program when they start to run, changes made will only take effect on the other stations when the program is re-started on that station.

---

## 17.5 Exclusive Processes

Some processes need to use the files exclusively. They are:

- Backing up data
- Restoring data
- Re-Indexing data
- Removing deleted records

For these processes to run, all the other users must either quit the program, or return to the program's main menu.

Some processes need to lock a file or a record. They are:

- Add
- Edit
- Delete
- Undelete
- Tag

---

## 17.6 Concurrent Processes

Some processes can operate concurrently with other processes, while others cannot.

Process 1	Concurrent Process
- Edit Employee A	- Edit Employee B
- View next, previous	- Edit/add
- Reports	- Edit/add
- Lists	- Edit/add
- Delete/undelete/tag Employee A	- Edit/add Employee B
- Backup	- None allowed
- Restore	- None allowed
- Re-index	- None allowed
- Remove deleted	- None allowed



## CHAPTER 18 - Interfacing with Other Softwares

### 18.1 Medical Director

---

“Medical Director” is a popular software for Medical Centre and Doctor’s surgeries. This software can be setup to look every minutes for a file containing one or several new patient entries and to add it to its database.

Medical Director Setup must be altered to setup the link and tell it where the file will be located. This information is available in Medical Director’s help files.

Workplace Audiometrics must also be setup in the following manner:

- Go to Setup > Preference and tick “Activate Link with MD”
- Click on the Path button and select where the file must be save

Once setup, Workplace Audiometrics will automatically place a file with the correct information where Medical Director expect the file. Every 60 seconds, Medical director will read that file and add the content to its database.

This function is specially useful if Clients are entered in Workplace Audiometrics before they are recorded in Medical Director. If they are already entered in Medical Director this function will not be useful as you would get duplicates.

Using the Audiogram in Reports see next paragraph.

Medical Director is a trademark of “HCN, Health Communication Network”

### 18.2 MS Word

---

You can include the graph displayed in the Employee data entry window in a Word or Excel document using the Windows Clipboard. When adding or editing a Test, click on the green button at the bottom of the graph to copy it to the Clipboard.

To paste the graph, move to the document, position the cursor where you want the graph inserted, and press **Ctrl V** or click on **Edit > Paste**.

In MS Word you can Right Click the image and alter its properties and behaviour.

### 18.3 MS Excel

---

The Employee and Test can be exported to Excel in List > Export to Excel. You will be prompted to enter a file name and a location. The file will be saved as a csv file. To open it in Excel, just select the file type .csv, you can then save it as an excel workbook if you wish. The same selection criteria as for other list apply.

### 18.4 SAP

---

This interface can be used for a wide range of purpose when you have the ability to read data and write a script to control the interface. It was designed initially to interface with SAP.

The interface work in the following way

Every time a records is added, edited or deleted, the Workplace Audiometrics checks for the presence of a

text file called SAP.txt, if it does not exist, the program creates it and stores the type of action (A=Add, E=Edit, D=Delete) followed by the data.

SAP (or any other Program) is meant to read the file and delete it from time to time. Workplace Audiometrics will then create a new one for the next and following actions, until it is deleted again. The deletion is required to keep the file from growing too big over time..

The interface needs to be activated in **Setup > Preference >Exporting data**, Tick the **Auto-save** checkbox.

The file is a text file, comma delimited and each field is enclosed in double quotes.

The format of the data is as follow:

- 1- Action (A=Add, E=Edit, D=Delete)
- 2- Surname
- 3- First Name
- 4- Employee number
- 5- Code
- 6- Sex
- 7- Date of Test
- 8- Visual examination Left
- 9- Left 500Hz
- 10- Left 1000Hz
- 11- Left 1500Hz
- 12- Left 2000Hz
- 13- Left 3000Hz
- 14- Left 4000Hz
- 15- Left 6000Hz
- 16- Left 8000Hz
- 17- Normal
- 18- Percentage loss Left
- 19- Visual examination Right
- 20- Right500Hz
- 21- Right1000Hz
- 22- Right1500Hz
- 23- Right2000Hz
- 24- Right3000Hz
- 25- Right4000Hz
- 26- Right6000Hz
- 27- Right8000Hz
- 28- Right
- 29- Percentage loss Right
- 30- Percentage loss Binaural
- 31- Occupation
- 32- Department
- 33- Shift



## CHAPTER 19 The To-Do List

To-Do's are meant to act as reminders of tasks to be performed at a certain time. They are not automatic reminders, you need to list your tasks for the day at the beginning of each day, or on a weekly basis. To-Do's are not exported to WA text files, nor imported.

The To-do list can be accessed from 2 different locations



### Within the Data Entry Screen

Click on the To-do Tab to display the To-Do's associated with an Employee.

This To-Do will also appear on the Global To-Do list available from the **Main** menu. The name of the Employee will show automatically on the Global To-Do but not in the **Data entry** To-Do. This To-Do is meant to record an action associated with an employee.



### From the Main Menu

This To-do is not connected with any employee, will not appear in the Data entry screen, but will appear in the Global To Do. It is meant to record a task not associated with a particular employee, like "Calibrate Audiometer" for example.

### 19.1 Adding a To-do

---



#### Employee To-Do

Go to Data Enter, Click on the To-Do Tab, Click on the **Add** button, select **date**, **time due** and a **description** of the task.



#### Global To Do

Go the **Main** menu, Click on the **To-Do** button, Click on the **Add** button, enter as for above.,

### 19.2 To-do Display

---

To-Do's display differently when they are due, overdue or completed.

To mark a To-do as **Completed**, just double-click on the **Done** check box within the grid

The To-dos overdue show in red, the completed ones show in default light colour.

Normally the Completed To-Do's do not show anymore, unless you specifically ask for it by ticking the **Show All** checkbox.

### 19.3 To-do Categories

---

The list of To-do categories can be added to in **Setup > To-Do Categories**. There is no limit to the number of categories

### 19.4 Listing the To-dos

---

To list To-Dos, go to **List > To-do** and make your selection of criteria.





## CHAPTER 20 The Security System


### 20.1 Enabling the Security System

To enable the security system go to Setup > Preference and tick the check box below.



The screenshot shows the 'Options' tab in the 'Setup > Preference' window. The 'Date Format' is set to 'D/M/Y' and 'Australia'. The 'Colour Scheme' is 'Classic (high res. blue)'. In the 'Passwords' section, the 'Enable Password Security' checkbox is checked, and a padlock icon is visible. A green arrow points to the checkbox.

Once this is enabled you will have an additional menu option in the Setup menu called **Security**



The screenshot shows the 'Data Access Management' window. It contains a table with the following data:

User name	Updated	Active	Employee Add	Edit	Delete	Assess	Test Add	Edit	Delete	Todo Add	Edit	Delete	Setup
ADMIN	13/08/2012	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
MASTER	13/08/2012	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SUPERUSER	25/07/2012	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
USER	13/08/2012	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Below the table are buttons: Add, Edit, Del/Undel, and Exit. A help box at the bottom contains the following information:

- Each Logon name has 6 rights which applies to all data entry windows.
- The Admin Logon name cannot be deleted or changed except for the password.
- The Superuser Logon is reserved and cannot be altered in anyway.
- A Menu Access Profile is connected to each logon name. It defines which menu items will be enabled for this logon.

In this screen you can add users and enable or disable any of their rights

Logon Details and Access Rights

Logon NameADMIN

Active☒

Password\*\*\*\*\*

Re-enter Password\*\*\*\*\*

Max length: 10 chr

Employee

☒ Add Employee

☒ Edit Employee

☒ Delete Employee

Test

☒ Add Test

☒ Edit Test (same day only)

☒ Edit Test (anytime)

☒ Delete Test

To-Do

☒ Add To-Do

☒ Edit To-Do

☒ Delete To-Do

Regardless of the Edit settings, Data can always be edited the same day it is entered

Only the password can be changed in this logon

Entered25/07/2012

Updated13/08/2012

Save

Cancel

General

☒ Utility 1

☒ Utility 2

☒ Setup Menu

☒ Preferences

☒ Security

Data

☒ Export Data

☒ Import Data

☒ Backup Data

☒ Restore Data

☒ Potential Liability

☒ Prevent printing from from Data for a specific employee

You can select from a number of access rights to give each use exactly what you want as an administrator.





## CHAPTER 21 Headset Correction

### 21.1 What is Headset Correction

Headset correction is a set of dB value which needs to be subtracted from each threshold, to correct for loss due to the shape of the cavity between the headphone insert and the ear of the person being tested. The amount to subtract is different for each frequency. It is often higher at low frequencies.


#### - When is it needed?

Normally the manufacturer will provide this information if it is required for their headset.

#### - Why is it needed?


When an audiometer is calibrated, the inserts inside the headphone are removed from the headphone and set against a very sensitive calibration microphone.

However when the insert is located inside the headphone, it may not sit right against the ear of the person being tested, like it did against the microphone during the calibration process. There may be a gap between the insert and the ear.



In Fig 1, the sound source is in the same position relative to the ear as it is to the microphone during the audiometer calibration. In this case there is no need for correction.

The Audiocup Headset belongs to this category.



In Fig 2, an additional distance between the sound source and the ear is present due to the headphone padding in gray.

The sound has further to travel to get to the ear (area in green).

The air gap induces a loss, as the sound waves have further to travel to get to the ear.

Depending on the manufacturer's instruction, there may be a need to introduce a correction. The Bilsom 847 and Peltor 2000 headsets belong to this category.

## 1.1 Enabling Headset Correction

---

The correction needs to be enabled in **Setup > Audiometer > Enable correction**

In most instance this is not necessary.

- Activating the correction for a specific audiometer

Once correction is enabled, it needs to be activated independently for each audiometer, as the correction values may be different for each and some may not require it.

- Go to **Setup > Audiometer**

- Edit the audiometer record

- Tick the **Activate Correction** checkbox

- Enter the correction factor for each frequency, zero means no correction required at that frequency

- Enter the name of the headset

- **Save**

### Correction data for common headsets

Peltor 2000 Bilsom 847		
.5k	10 dB	10 dB
1k	5 dB	5 dB
1.5k	5 dB	5 dB
2k	5 dB	5 dB
3k	5 dB	10 dB
4k	5 dB	5 dB
6k	10 dB	5 dB
8k	0 dB	0 dB

Values are entered in 5dB steps

### When is correction applied

Once correction is enabled, and activated for a specific audiometer a new button to apply the correction to a test result will appear under the heading "**Headset**" in the test screen.

#### - When Capturing a test

When capturing a test, and correction is enabled, and activated the test will be corrected automatically, even if you re-capture over a test. Following this, the Correction button will be disabled and a green led light will be present near the correction button. If you edit that test, the correction button will no longer be visible.

However if you **Re-Capture** the test, it will be corrected also.

If the Australian standard is selected, the retest is corrected in the same manner.

#### - When Entering a test manually

After you have entered a test manually, you need apply the correction by clicking on the Correction button. You can only do this once as the button will then be disabled. The only way to reverse this is by manually entering the results again.

## APPENDIX A - Layout of Export/Import file

The import function expect an ASCII text file with the following fields. Character fields are enclosed in double quotes, Numeric fields are not. Each field is separated by a comma, each line is separated by a carriage return-Line feed

Ex "V3.0,"BLOGGS","fred",28,"12/12/2000"

In this example, 28 is a numeric field, the other fields are characters

1	VERSION	Character	4	"V3.0", or any dummy field (see Note below)
2	SURNAME	Character	18	
3	NAME	Character	18	
4	SEX	Character	1	
5	FILENO	Character	8	Employee number
6	CODE	Character	10	
7	STATUS	Character	1	Employee status
8	STREET	Character	28	
9	CITY	Character	20	
10	STATE	Character	3	
11	ZIP	Character	4	
12	TEL	Character	13	
13	DOB	Character	8	
14	COM_DATE	Character	8	Commencement date
15	ENTERED	Character	8	
16	RECALL	Character	8	Re-test date
17	RETRAIN	Character	8	Re-train date
18	TERMINATED	Character	8	
19	RECALLYES	Character	1	Y/N
20	RETRAINYES	Character	1	Y/N
21	SHIFT	Character	5	
22	CLASS	Character	27	Occupation
23	DEPARTMENT	Character	27	
24	WORKGROUP	Character	12	
25	NOISECLASS	Numeric	1	Workgroup noise level
26	PROTCLASS	Numeric	1	Hearing Protector class
27	COMPANY	Character	30	
28	ADDRESS	Character	30	
29	COMMENT1	Character	51	
30	PREEMP1	Character	36	Previous employer 1
31	PREEMP2	Character	36	
32	PREEMP3	Character	36	
33	MILITARY	Character	1	Y/N
34	TIMESERVED	Character	2	Years
35	EXPLOSIVES	Character	1	Y/N
36	TROUBLE	Character	1	Y/N
37	HOBBIES	Character	1	Y/N
38	POWERTOOLS	Character	1	Y/N
39	DIZZINESS	Character	1	Y/N
40	RINGING	Character	1	Y/N
41	COLDS	Character	1	Y/N
42	OTHER	Character	1	Y/N (discharges)
43	REFERRED	Character	1	Y/N
44	ANTIBIOTIC	Character	1	Y/N
45	MUMPS	Character	1	Y/N
46	HISTCHANGE	Character	1	Y/N
47	COMMENT2	Character	130	
48	NBTEST	Numeric	2	Number of tests
49	TAG	Numeric	1	
50	PERS	Numeric	1	
51	SPARE1	Character	1	
52	SPARE2	Character	1	
53	AGE	Character	2	
54	DATE_TEST	Character	8	
55	TRECALL	Character	8	re-test date
56	TRETRAIN	Character	8	
57	LEN_EMPLOY	Numeric	2	
58	TCLASS	Character	27	occupation

59	TDEPT	Character	27	department
60	QUIET	Character	3	period of quiet before test
61	PROTWORN	Character	1	protection worn
62	PROTUWORN	Character	1	protection worn before test
63	PROTOTYPE	Character	16	type of protector
64	TPROTCCLASS	Numeric	1	protector class
65	TRAINED	Character	1	Y/N
66	TESTTYPE	Character	1	
67	MAKE1	Character	15	Test instrument
68	CAL1	Character	10	Calibration date
69	SN1	Character	7	Serial no
70	AUDIOLOG	Character	18	Testing Officer
71	APPRO	Character	8	Approval number
72	VISUAL_L	Character	8	Visula examination
73	VISUAL_R	Character	8	
74	AIRNORM_L	Character	1	
75	AIRNORM_R	Character	1	
76	ABNORMAL	Character	1	Abnormal test
77	PER_LEFT	Character	5	% loss
78	PER_RIGHT	Character	5	
79	BINAURAL	Character	5	
80	TCOMMENT11	Character	250	Test comment part 1
81	TCOMMENT12	Character	110	Test comment part 2
82	AIRC_TEST	Character	32	Air test
83	RE_TEST	Character	32	Re test
84	BONE_TEST	Character	32	Bone test
85	LMEP	Character	4	Tymp
86	RMEP	Character	4	
87	LPV	Character	4	
88	RPV	Character	4	
89	LCOMP	Character	5	
90	RCOMP	Character	5	
91	LGRAD	Character	3	
92	RGRAD	Character	3	
93	FCV	Character	4	Spirometry
94	FEV1	Character	4	
95	RATIO	Character	3	
96	FEF	Character	4	
97	PEF	Character	5	
98	TEARACHE	Character	1	Y/N
99	TEXPLOSIVE	Character	1	Y/N
100	THOBBIES	Character	1	Y/N
101	TPOWERTOOL	Character	1	Y/N
102	TDIZZINESS	Character	1	Y/N
103	TRINGING	Character	1	Y/N
104	TCOLDS	Character	1	Y/N
105	TREFERRED	Character	1	Y/N
106	TOTHER	Character	1	Y/N (discharges)
107	TSPARE1	Character	1	
108	TSPARE2	Character	1	
109	SPARE3	Character	1	
110	SPARE4	Character	1	
111	C_ADD2	Character	30	
112	C_CITY	Character	20	
113	C_STATE	Character	8	
114	C_ZIP	Character	8	
115	C_TEL	Character	13	
116	C_FAX	Character	13	
117	C_EMAIL	Character	60	
118	C_CONTACT	Character	30	
119	C_EMPLOYER	Character	30	

The Employee Status options are as follow:

P- Permanent noise exposed	N-Permanent quiet environment	W- Contractor	A- Applicant
H- Permanent high risk	T- Terminated	O- Trainee	D- Deceased
U- Other			C- Casual

**Note :** If the first field does not contain "V3.0" or "V4.0" the software will issue a warning that the type of data is unknown. If you loaded this field with dummy data, this is fine, you can ignore the warning, the function will proceed as normal. When exporting from a mainframe for example, it is easier to put a dummy field in this position.

## APPENDIX B - Layout of Update file

The data layout is the same as for the Import/Export file, but only the first 28 fields are required. If the rest is present, it is ignored.

If the Department or the Company in the update file do not exist in the program's list of Departments or Companies, the new Department or Company is added to that list. If the Department or Company already exists, the number is updated automatically.

The Date Entered is replaced with the date of the Update

Only the fields with a ✓ are updated, however all the fields up to field 28 (as listed in the import file layout) are required in their proper places. A field can be empty if no change is required ("").

Only the fields in the update file containing data are used for updating the database. If a field is empty in the update file and the corresponding field in the database contains data, this data will not be altered.

1	VERSION	Character	4	"V3.0", or dummy field, (see note below)
2	SURNAME	Character	18	
3	NAME	Character	18	
4	SEX	Character	1	
5	FILENO	Character	8	Employee number
✓ 6	CODE	Character	10	
✓ 7	STATUS	Character	1	Employee status
✓ 8	STREET	Character	28	
✓ 9	CITY	Character	20	
✓ 10	STATE	Character	3	
✓ 11	ZIP	Character	4	
✓ 12	TEL	Character	13	
13	DOB	Date	8	
14	COM_DATE	Date	8	Commencement date
15	ENTERED	Date	8	
✓ 16	RECALL	Date	8	Re-test date
✓ 17	RETRAIN	Date	8	Re-train date
✓ 18	TERMINATED	Date	8	
✓ 19	RECALLYES	Character	1	Y/N
✓ 20	RETRAINYES	Character	1	Y/N
✓ 21	SHIFT	Character	5	
✓ 22	CLASS	Character	27	Occupation
✓ 23	DEPARTMENT	Character	27	
✓ 24	WORKGROUP	Character	12	
✓ 25	NOISECLASS	Numeric	1	Workgroup noise level
✓ 26	PROTCLASS	Numeric	1	Hearing Protector class
✓ 27	COMPANY	Character	30	
✓ 28	ADDRESS	Character	30	

To make a sample file, tag a couple of Employees then go to LIST, select EXPORT DATA, select for ALL TAGGED. This will give you a small file to play with, .

**Note:** If the first field does not contain "V3.0" the software will issue a warning that the type of data is of unknown type. This is to protect from importing data from the DOS version of this software which does not have the same fields.

If you know this field is dummy data, ignore the warning, the process will proceed as normal.

When exporting from a mainframe for example, it is easier to put a dummy field in this position. If this field was omitted you will damage the datafile and fill it with garbage.

**Precautions:**

It is recommended to do a **backup** before this operation if in doubt about the validity of the data to allow for rollback (using the **Utility > Restore** function).



## INDEX

Live Spell 74  
Windows 95 5  
600 x 800 7  
Actions Allowed 45  
AD27 60  
Add 29, 81, 82  
Allow Re-Test 15  
Allow Re-Train 15  
AS53 66  
AS54 66  
Audiogram 32, 35  
Australian Standard 11  
Australian Standard 1269 50  
Backup 30, 48, 82  
Baseline test 45  
Batch 32, 35  
Calendar 26  
Card 35  
Client letter 70  
Combo Box 25  
Comment 55  
Company 10  
Control Panel 41  
Copyright 5  
data source 43  
Date Field 26  
Delete 35, 47, 82  
Deleted 31  
Delimited 42  
Department 10  
Desktop area 7  
Dictionaries 74  
Earscan 62  
Edit 29, 35, 82  
Edit window 25  
Envelope 73  
Excel 44  
Field names 57  
Fields 42  
File 29, 31, 47  
Files 29, 82  
Find 35, 47  
Form 27 20, 32, 35  
Format 42  
format(Rich Text 20  
Graph Colors 54  
Grid 26  
Hearing Protection 46  
index 47  
Indexing 82  
Labels 44  
Letters 42  
Live Spell 74  
MA728 64  
MA800 62  
Mailmerge 42, 57  
Medical History 15  
Merge fields 44  
MICROLAB 62  
Microsoft Word 42  
MS Word 20, 42  
NAL 50  
network 81, 83  
New tables (1988) 53  
Normals 77  
Occupation 10  
Old table (1973) 53  
OSCILLA SM910 59, 61, 67  
OSCILLA SM930 59, 61, 67  
OSHA 36, 37  
Page Setup 41  
Password 49, 55  
Percentage loss 11, 15, 18, 19, 76, 78, 79  
Potential Liability 38  
Pre-View Reports 19  
Presbycusis 76  
Purge deleted 47  
Re-index 47  
Registration 5  
Replaceable field 70  
Report Card 19  
SCREENW.INI 48  
Shortcut 7, 9  
Shortcut 7  
SLC80 53  
Speed button 25  
Speed Buttons 11  
Tabbed Book 25  
tag 47  
Test Instrument 16  
Threshold Normals 54  
Tonal noise 53  
Updates 5  
Windows 98 7  
WinZip 31, 48  
Word Perfect 42  
Word processor 42  
WordStar 42  
Work number 29  
Workgroup 10  
Workstation 9  
ZIP. 30

-